

Pointe Scientific Inc.
Stat Fax
Application Settings
Rev: 4-02

Acid Phos

Rate By Factor
Wavelengths: 405 600 nm
Lag Time: 120 Secs
Read Time: 60 Secs
Factor: 1628
Norms: 0.0 to 9.0
Linear: 0.0 to 35.0
Sample Volume: 50 ul
Reagent Volume: 1 ml

Albumin

Standard
Wavelengths: 600 490 nm
Standard #1: User Defined
Norms: 3.5 to 5.3
Linear: 0.0 to 6.0
Sample Volume: 10 ul
Reagent Volume: 1 ml

Alk. Phos

Rate By Factor
Wavelengths: 405 600 nm
Lag Time: 30 Secs
Read Time: 30 Secs
Factor: 2187
Norms: 35 to 123
Linear: 0.0 to 800
Sample Volume: 25 ul
Reagent Volume: 1 ml

ALT

Rate By Factor
Wavelengths: 340 405 nm
Lag Time: 60 Secs
Read Time: 60 Secs
Factor: -3376
Norms: 0.0 to 38.0
Linear: 0.0 to 700
Sample Volume: 50 ul
Reagent Volume: 1 ml

Amylase

Rate By Factor
Wavelengths: 405 600 nm
Lag Time: 60 Secs
Read Time: 60 Secs
Factor: 3178
Norms: 25 to 125
Linear: 0.0 to 2000
Sample Volume: 25 ul
Reagent Volume: 1 ml

AST

Rate By Factor
Wavelengths: 340 405 nm
Lag Time: 60 Secs
Read Time: 60 Secs
Factor: -3376
Norms: 0.0 to 40.0
Linear: 0.0 to 700.0
Sample Volume: 50 ul
Reagent Volume: 1ml

Total Bilirubin

Standard
Wavelengths: 550 600 nm
Standard: User Defined
Norms: 0.2 to 1.2
Linear: 0.0 to 20.0
Sample Volume: 50 ul
Reagent Volume: 1 ml

Direct Bilirubin

Factor
Wavelengths: 550 600 nm
Factor: 45.08
Norms: 0.0 to 0.5
Linear: 0.0 to 20.0
Sample Volume: 100 ul
Reagent Volume: 1 ml

BUN

Rate By Standard
Wavelengths: 340 405 nm
Standard: User Defined
Lag Time: 30 Secs
Read Time: 60 Secs
Norms: 7.0 to 18.0
Linear: 0.0 to 80.0
Sample Volume: 10 ul
Reagent Volume: 1 ml

Calcium

Factor
Wavelengths: 550 600 nm
Standard: User Defined
Norms: 8.5 to 10.4
Linear: 0.0 to 20.0
Sample Volume: 25 ul
Reagent Volume: 1 ml

Chloride

Standard
Wavelengths: 490 600 nm
Standard: User Defined
Norms: 98.0 to 106.0
Linear: 70.0 to 130.0
Sample Volume: 10 ul
Reagent Volume: 1 ml

Cholesterol

Standard
Wavelengths: 490 600nm
Standard: User Defined
Norms: 120 to 240
Linear: 0 to 500
Sample Volume: 10 ul
Reagent Volume: 1 ml

CPK

Rate By Factor
Wavelengths: 340 405 nm
Lag Time: 120 Secs
Read Time: 30 Secs
Factor: 6592
Norms: 25 to 192
Linear: 0 to 1500
Sample Volume: 25 ul
Reagent Volume: 1 ml

CO2

Standard
Wavelengths: 340 405 nm
Standard: User Defined
Norms: 22.0 to 29.0
Linear: 0.0 to 40.0
Sample Volume: 10
Reagent Volume: 1 ml

Creatinine

Rate By Standard
Wavelengths: 490 600 nm
Standard: User Defined
Lag Time: 20 Secs
Read Time: 60 Secs
Norms: 0.7 to 1.7
Linear: 0.5 to 7.0
Sample Volume: 50 ul
Reagent Volume: 1 ml

Hemoglobin

Standard
Wavelengths: 550 600 nm
Standard: User Defined
Norms: 11.0 to 18.0
Linear: 0.0 to 20.0
Sample Volume: 5 ul
Reagent Volume: 1 ml

GGT

Rate By Factor
Wavelengths: 405 600 nm
Factor: 4316
Lag Time: 30 Secs
Read Time: 30 Secs
Norms: 8.0 to 54.0
Linear: 0.0 to 700.0
Sample Volume: 25 ul
Reagent Volume: 1 ml

Glucose Hex

Standard
Wavelengths: 340 405 nm
Standard: User Defined
Norms: 65 to 110
Linear: 0.0 to 500
Sample Volume: 10 ul
Reagent Volume: 1 ml

Glucose Ox

Standard
Wavelengths: 490 600 nm
Standard: User Defined
Norms: 65 to 110
Linear: 0.0 to 500
Sample Volume: 10 ul
Reagent Volume: 1 ml

HDL Cholesterol

Standard
Wavelengths: 490 600nm
Standard: User Defined
Norms: 30 to 75
Linear: 0 to 200
Sample Volume: 25 ul
Reagent Volume: 1 ml

Iron

Rate By Standard
Wavelengths: 550 600 nm
Lag Time: 10 Secs
Read Time: 60 Secs
Standard: User Defined
Norms: 60 to 150
Linear: 0 to 500
Sample Volume: 200 ul
Reagent Volume: 1 ml buffer, 20 ul color

LD-L

Rate By Factor
Wavelengths: 340 405 nm
Factor: 6592
Lag Time: 60 Secs
Read Time: 30 Secs
Norms: 80 to 227
Linear: 0 to 800
Sample Volumes: 25 ul
Reagent Volumes: 1 ml

Lipase

Rate By Standard (Batch Mode)
Wavelengths: 550 600 nm
Standard: User Defined
Lag Time: 180 Secs
Read Time: 280 Secs
Norms: 0 to 62
Linear: 0 to 600
Sample Volume: 10 ul
Reagent Volume: 1 ml

Magnesium

Standard
Wavelengths: 550 600 nm
Standard: User Defined
Norms: 1.3 to 2.5
Linear: 0.0 to 4.0
Sample Volume: 10 ul
Reagent Volume: 1 ml

Phosphorus

Standard
Wavelengths: 340 405 nm
Standard: User Defined
Norms: 2.5 to 4.8
Linear: 0.0 to 12.0
Sample Volume: 10 ul
Reagent Volume: 1 ml

Potassium

Standard
Wavelengths: 600 nm
Standard: User Defined
Norms: 3.4 to 5.3
Linear: 2.0 to 6.5
Sample Volume: 25 ul
Reagent Volume: 1 ml

Protein

Standard
Wavelengths: 550 600 nm
Standard: User Defined
Norms: 6.2 to 8.5
Linear: 1.0 to 15
Sample Volume: 25 ul
Reagent Volume: 1 ml

Triglyceride

Standard
Wavelengths: 550 405 nm
Standard: User Defined
Norms: 36 to 165
Linear: 0 to 1000
Sample Volume: 10 ul
Reagent Volume: 1 ml

Uric Acid

Standard
Wavelengths: 490 600 nm
Standard: User Defined
Norms: 2.5 to 7.7
Linear: 0.0 20.0
Sample Volume: 25 ul
Reagent Volume: 1 ml

UIBC

Rate By Standard
Wavelengths: 550 600 nm
Standard: 500
Lag Time: 10 Secs
Read Time: 60 Secs
Norms: User Defined
Linear: 0.0 to 500.0
Sample Volumes: 125 ul
Reagent Volumes: 1 ml buffer, 25 ul color

Micro Protein

Standard
Wavelengths: 600 nm
Standard: User Defined
Norms: 0.0 to 10.0
Linear: 0 to 250.0
Sample Volume: 20 ul
Reagent Volume: 1 ml

CK- MB

Rate By Factor
Wavelengths: 340 405 nm
Lag Time: 300 Secs
Read Time: 60 Secs
Factor: 6752
Norms: 0 to 22
Linear: 0 to 2000
Sample Volume: 50 ul
Reagent Volume: 1 ml

Sodium

Standard
Wavelengths: 550 nm
Standard: User Defined
Norms: 135 to 155
Linear: 0 to 200
Refer to Package insert for
sample and reagent volumes