

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ACID PHOSPHATASE
Catalog # : A7503

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	U/L

Analysis:

P SAMPLE:	20 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	280.0 S
NUMBER:	29
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI:40.0 UL
NORMAL RANGE:	LO: 0.0 HI:9.0 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 29
REAC. LIMIT:	0.3000
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	969.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ALKALINE PHOSPHATASE
Catalog # : A7505,A7516

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	150 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:800 UL
NORMAL RANGE:	LO: 35 HI:123 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.7500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	2717
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:1.000
BLANK RANGE:	LO: -0.005 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ALT
Catalog # : A7525, A7526

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	16 UL
DILUENT:	30 UL
REAGENT:	145 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	29
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	DECREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 UL
NORMAL RANGE:	LO: 0 HI:38 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 29
REAC. LIMIT:	0.2500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH RUN
CALIBRATOR:	
FACTOR:	2512.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
 Test: AST
 Catalog # : A7560,A7561

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	16 UL
DILUENT:	30 UL
REAGENT:	145 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	15
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	DECREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0	HI:500 UL
NORMAL RANGE:	LO: 0	HI:40 UL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 15
REAC. LIMIT:	0.2500	
POINT:	M1	

Calibration:

CALIB. INTERVAL:	EACH RUN	
CALIBRATOR:		
FACTOR:	2512.00	
POS:		
REPLICATE:		
REAGENT RANGE:	LO: NO	HI:NO
BLANK RANGE:	LO: NO	HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: BUN
Catalog # : B7550,B7552

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	60 UL
REAGENT:	300 UL
AO READINGS	
FIRST:	20.0 S
NUMBER:	10
INTERVAL:	10 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	DECREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI:80.0 MG/DL
NORMAL RANGE:	LO: 7.0 HI:18.0 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 10
REAC. LIMIT:	NO
POINT:	

Calibration:

CALIB. INTERVAL:	
CALIBRATOR:	20.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 1.000 HI:3.000
BLANK RANGE:	LO: -0.100 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CALCIUM (ARSENAZO)
Catalog # : C7529

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	630 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	2 UL
DILUENT:	50 UL
REAGENT:	150 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	40 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:15.00 MG/DL
NORMAL RANGE:	LO: 8.5 HI:10.4 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	10.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: -0.100 HI:2.000
BLANK RANGE:	LO: 0.100 HI:2

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CHOLESTEROL
Catalog # : C7509, C7510

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	500 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3.0 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 MG/DL
NORMAL RANGE:	LO: 0 HI:240
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	200 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CPK
Catalog # : C7512, C7522

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1500 UL
NORMAL RANGE:	LO: 25 HI:192 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.3200
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	8039.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.080 HI:0.380
BLANK RANGE:	LO: -0.010 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: CREATININE

Catalog # : C7539

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	520 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	20 UL
DILUENT:	40 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	35.0 S
NUMBER:	10
INTERVAL:	10 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0.0	HI:25.0 MG/DL
NORMAL RANGE:	LO: 0.7	HI:1.5 MG/DL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 10
REAC. LIMIT:		
POINT:		NO

Calibration:

CALIB. INTERVAL:	EACH TEST	
CALIBRATOR:	5.0 MG/DL	
FACTOR:		
POS:	3	
REPLICATE:	SINGLE	
REAGENT RANGE:	LO: 0.025	HI:0.200
BLANK RANGE:	LO: -0.003	HI:0.003

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: GGT

Catalog # : G7571

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	10 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0	HI:600 UL
NORMAL RANGE:	LO: 9	HI:54 UL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 16
REAC. LIMIT:	0.2500	
POINT:	M1	

Calibration:

CALIB. INTERVAL:	EACH TEST	
CALIBRATOR:		
FACTOR:	2525.00	
POS:		
REPLICATE:		
REAGENT RANGE:	LO: 0.100	HI:1.000
BLANK RANGE:	LO: -0.001	HI:0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: GGT (SOLUBLE)
Catalog # : G7570

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	10 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1000 UL
NORMAL RANGE:	LO: 9 HI:54 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.2500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	2631.0
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:1.000
BLANK RANGE:	LO: -0.001 HI:0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: GLUCOSE HEXOKINASE
Catalog # : G7518, G7517

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	20 UL
REAGENT:	230 UL
AO READINGS	
FIRST:	5.0 S
NUMBER:	2
INTERVAL:	120 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 MG/DL
NORMAL RANGE:	LO: 65 HI:110 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	100 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.0000 HI:0.250
BLANK RANGE:	LO: -0.015 HI:0.005

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: LDH
Catalog # : L7535, L7572

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:800 UL
NORMAL RANGE:	LO: 80 HI:285 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 5 LAST: 16
REAC. LIMIT:	0.1700
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	8039.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:0.600
BLANK RANGE:	LO: -0.005 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: Lp(a)
Catalog # : L7597

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-A0
CALIBRATION MODE:	STD Non Linear
REAGENT BLANK:	RGT/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	11UL
DILUENT:	30 UL
REAGENT:	210 UL
I1 Incubation	240 S
M1	220 S
SR1 Start Reagent 1	35 S
Diluent	20 S
AO READINGS	
FIRST:	5 S
NUMBER:	6 S
INTERVAL:	240 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	OFF	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: NO	HI:NO
NORMAL RANGE:	LO: 0	HI:30
CALC. STEPS:	1	
CALC. STEP A:	Endpoint	
READINGS:	FIRST: M1	LAST: 6
REAC. LIMIT:		
POINT:		

Calibration:

CALIB. INTERVAL: On Request	Replicate : 2 Duplicate
Deviation %: 10	Calculation Model : Linear Interpolation

Calibration Correction : NO

Antigen Excess : NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: TRIGLYCERIDES (GPO)
Catalog # : T7531, T7532

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	540 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1000 MG/DL
NORMAL RANGE:	LO: 36 HI:165 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	300 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.000 HI:0.200
BLANK RANGE:	LO: -0.020 HI:0.000

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: URIC ACID
Catalog # : U7580, U7581

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	520 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	5 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0.0	HI:25.0 MG/DL
NORMAL RANGE:	LO: 1.5	HI:7.0 MG/DL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 2
REAC. LIMIT:		
POINT:		

Calibration:

CALIB. INTERVAL:	EACH TEST	
CALIBRATOR:	5.0 MG/DL	
FACTOR:		
POS:	3	
REPLICATE:	SINGLE	
REAGENT RANGE:	LO: 0.000	HI:0.410
BLANK RANGE:	LO: -0.020	HI:0.050

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: ALBUMIN

Catalog # : A7502

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	630 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	G/DL

Analysis:

R1 REAGENT 1:	370 UL
I1 INCUBATION:	15 S
M1:	10 S
S SAMPLE:	2 UL
DILUENT:	28 UL

AO READINGS

FIRST:	0.5 S
NUMBER:	2
INTERVAL:	60 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.5G/DL HI: 6.0 G/DL
NORMAL RANGE:	LO: 3.5 G/DL HI: 5.3 G/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR CAL-1:	4.0 G/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.600 HI: 0.800
BLANK RANGE:	LO: 0.060 HI: 1.100

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: PHOSPHORUS
Catalog # : P7516

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

R1 REAGENT 1:	200 UL
I1 INCUBATION:	5 S
M1:	4 S
S SAMPLE:	8 UL
DILUENT:	30 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	90 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 12.0
NORMAL RANGE:	LO: 2.5 HI: 4.8 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR CAL-1:	5.0 MG/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: NO HI: NO
BLANK RANGE:	LO: NO HI: NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: TOTAL PROTEIN

Catalog # : T7528

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	540 NM
TEMPERATURE:	37C
DECIMAL POSITION:	2
UNIT:	G/DL

Analysis:

R1 REAGENT 1:	200
I1 INCUBATION:	5 S
M1:	4 S
S SAMPLE:	4 UL
DILUENT:	30 UL

AO READINGS

FIRST:	0.5 S
NUMBER:	2
INTERVAL:	115 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 1 HI: 15.0 G/DL
NORMAL RANGE:	LO: 6.2 HI: 8.5 G/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH RUN
CALIBRATOR CAL-1:	8.0 G/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.050 HI: 0.200
BLANK RANGE:	LO: -0.001 HI: 0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test:

Catalog # :

General:

MEASUREMENT MODE:
REACTION MODE:
CALIBRATION MODE:
REAGENT BLANK:
WAVELENGTH:
TEMPERATURE:
DECIMAL POSITION:
UNIT:

Analysis:

R1 REAGENT 1:
I1 INCUBATION:
M1:
S SAMPLE:
DILUENT:

AO READINGS
FIRST:
NUMBER:
INTERVAL:

Calculation:

CONVERS. FACTOR:
OFFS:
REAC. DIRECTION:
CHECK:
SAMPLE LIMIT:
TEST RANGE: LO: HI:
NORMAL RANGE: LO: HI:
CALC. STEPS:
CALC. STEP A:
READINGS: FIRST: LAST:
REAC. LIMIT:
POINT:

Calibration:

CALIB. INTERVAL:
CALIBRATOR CAL-1:
POS:
REPLICATE:
REAGENT RANGE: LO: HI:
BLANK RANGE: LO: HI:

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: DIRECT BILIRUBIN
Catalog # : B7538

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	546 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	180 S
I1 INCUBATION:	185 S
P SAMPLE:	20 UL
DILUENT:	30 UL
SR 1 START R1:	20 UL
REAGENT:	100 UL
DILUENT:	30 UL
AO READINGS	
FIRST:	1.0 S
NUMBER:	2
INTERVAL:	300 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 0.0 HI: 0.5 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	269.0
POS:	
REPLICATE:	
REAGENT RANGE:	LO: -0.0150 HI: 0.0050
BLANK RANGE:	LO: -0.0050 HI: 0.015

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: TOTAL BILIRUBIN
Catalog # : B7576

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	546 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	10 S
I1 INCUBATION:	15 S
P SAMPLE:	20 UL
DILUENT:	30 UL
SR 1 START R1:	20 UL
REAGENT:	100 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	1.0 S
NUMBER:	2
INTERVAL:	180 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 0.2 HI: 1.2 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	269.800
POS:	
REPLICATE:	
REAGENT RANGE:	LO: -0.0150 HI: 0.0050
BLANK RANGE:	LO: -0.0050 HI: 0.015

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CALCIUM
Catalog # : C7503, C7508

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	573 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	35 S
I1 INCUBATION:	40 S
P SAMPLE:	10 UL
DILUENT:	50 UL
SR 1 START R1:	50 UL
REAGENT:	150 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	40 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 8.5 HI: 10.4 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	10.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: -0.100 HI: 0.500
BLANK RANGE:	LO: 0.100 HI: 0.400

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CHLORIDE
Catalog # : C7501

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	480 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MEQ/L

Analysis:

M1:	120 S
I1 INCUBATION:	125 S
P SAMPLE:	5 UL
DILUENT:	25 UL
SR 1 START R1:	50 UL
REAGENT:	150 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	1
INTERVAL:	5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 80 HI: 120 MEQ/L
NORMAL RANGE:	LO: 98 HI: 106 MEQ/L
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 1

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	100 MEQ/L
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.0000 HI: 0.3000
BLANK RANGE:	LO: -0.010 HI: 0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
 Test: ACID PHOSPHATASE
 Catalog # : A7503

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	U/L

Analysis:

P SAMPLE:	20 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	280.0 S
NUMBER:	29
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI:40.0 UL
NORMAL RANGE:	LO: 0.0 HI:9.0 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 29
REAC. LIMIT:	0.3000
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	969.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ALKALINE PHOSPHATASE
Catalog # : A7505,A7516

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	150 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:800 UL
NORMAL RANGE:	LO: 35 HI:123 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.7500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	2717
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:1.000
BLANK RANGE:	LO: -0.005 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ALT
Catalog # : A7525, A7526

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	16 UL
DILUENT:	30 UL
REAGENT:	145 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	29
INTERVAL:	5 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	DECREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 UL
NORMAL RANGE:	LO: 0 HI:38 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 29
REAC. LIMIT:	0.2500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH RUN
CALIBRATOR:	
FACTOR:	2512.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
 Test: AST
 Catalog # : A7560,A7561

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	16 UL
DILUENT:	30 UL
REAGENT:	145 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	15
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	DECREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0	HI:500 UL
NORMAL RANGE:	LO: 0	HI:40 UL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 15
REAC. LIMIT:	0.2500	
POINT:	M1	

Calibration:

CALIB. INTERVAL:	EACH RUN	
CALIBRATOR:		
FACTOR:	2512.00	
POS:		
REPLICATE:		
REAGENT RANGE:	LO: NO	HI:NO
BLANK RANGE:	LO: NO	HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: BUN
Catalog # : B7550,B7552

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	60 UL
REAGENT:	300 UL
AO READINGS	
FIRST:	20.0 S
NUMBER:	10
INTERVAL:	10 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	DECREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI:80.0 MG/DL
NORMAL RANGE:	LO: 7.0 HI:18.0 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 10
REAC. LIMIT:	NO
POINT:	

Calibration:

CALIB. INTERVAL:	
CALIBRATOR:	20.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 1.000 HI:3.000
BLANK RANGE:	LO: -0.100 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CALCIUM (ARSENAZO)
Catalog # : C7529

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	630 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	2 UL
DILUENT:	50 UL
REAGENT:	150 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	40 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:15.00 MG/DL
NORMAL RANGE:	LO: 8.5 HI:10.4 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	10.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: -0.100 HI:2.000
BLANK RANGE:	LO: 0.100 HI:2

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CHOLESTEROL
Catalog # : C7509, C7510

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	500 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3.0 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 MG/DL
NORMAL RANGE:	LO: 0 HI:240
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	200 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: NO HI:NO
BLANK RANGE:	LO: NO HI:NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CPK
Catalog # : C7512, C7522

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1500 UL
NORMAL RANGE:	LO: 25 HI:192 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.3200
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	8039.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.080 HI:0.380
BLANK RANGE:	LO: -0.010 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: CREATININE

Catalog # : C7539

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	520 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	20 UL
DILUENT:	40 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	35.0 S
NUMBER:	10
INTERVAL:	10 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI:25.0 MG/DL
NORMAL RANGE:	LO: 0.7 HI:1.5 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 10
REAC. LIMIT:	
POINT:	NO

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	5.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.025 HI:0.200
BLANK RANGE:	LO: -0.003 HI:0.003

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: GGT

Catalog # : G7571

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	10 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0	HI:600 UL
NORMAL RANGE:	LO: 9	HI:54 UL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 16
REAC. LIMIT:		0.2500
POINT:		M1

Calibration:

CALIB. INTERVAL:	EACH TEST	
CALIBRATOR:		
FACTOR:	2525.00	
POS:		
REPLICATE:		
REAGENT RANGE:	LO: 0.100	HI:1.000
BLANK RANGE:	LO: -0.001	HI:0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: GGT (SOLUBLE)
Catalog # : G7570

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	405 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	10 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1000 UL
NORMAL RANGE:	LO: 9 HI:54 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 16
REAC. LIMIT:	0.2500
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	2631.0
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:1.000
BLANK RANGE:	LO: -0.001 HI:0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: GLUCOSE HEXOKINASE
Catalog # : G7518, G7517

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	20 UL
REAGENT:	230 UL
AO READINGS	
FIRST:	5.0 S
NUMBER:	2
INTERVAL:	120 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:500 MG/DL
NORMAL RANGE:	LO: 65 HI:110 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	100 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.0000 HI:0.250
BLANK RANGE:	LO: -0.015 HI:0.005

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: LDH
Catalog # : L7535, L7572

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A1
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	U/L

Analysis:

P SAMPLE:	5 UL
DILUENT:	30 UL
REAGENT:	200 UL
AO READINGS	
FIRST:	50.0 S
NUMBER:	16
INTERVAL:	10 S
M1:	0.5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:800 UL
NORMAL RANGE:	LO: 80 HI:285 UL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 5 LAST: 16
REAC. LIMIT:	0.1700
POINT:	M1

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	8039.00
POS:	
REPLICATE:	
REAGENT RANGE:	LO: 0.100 HI:0.600
BLANK RANGE:	LO: -0.005 HI:0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: TRIGLYCERIDES (GPO)
Catalog # : T7531, T7532

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	540 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	3 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0 HI:1000 MG/DL
NORMAL RANGE:	LO: 36 HI:165 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: 1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	300 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.000 HI:0.200
BLANK RANGE:	LO: -0.020 HI:0.000

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: URIC ACID
Catalog # : U7580, U7581

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-T-A0
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	520 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

P SAMPLE:	5 UL
DILUENT:	35 UL
REAGENT:	165 UL
AO READINGS	
FIRST:	3 S
NUMBER:	2
INTERVAL:	300 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	ON	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: 0.0	HI:25.0 MG/DL
NORMAL RANGE:	LO: 1.5	HI:7.0 MG/DL
CALC. STEPS:		
CALC. STEP A:		
READINGS:	FIRST: 1	LAST: 2
REAC. LIMIT:		
POINT:		

Calibration:

CALIB. INTERVAL:	EACH TEST	
CALIBRATOR:	5.0 MG/DL	
FACTOR:		
POS:	3	
REPLICATE:	SINGLE	
REAGENT RANGE:	LO: 0.000	HI:0.410
BLANK RANGE:	LO: -0.020	HI:0.050

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: ALBUMIN
Catalog # : A7502

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	630 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	G/DL

Analysis:

R1 REAGENT 1:	370 UL
I1 INCUBATION:	15 S
M1:	10 S
S SAMPLE:	2 UL
DILUENT:	28 UL

AO READINGS

FIRST:	0.5 S
NUMBER:	2
INTERVAL:	60 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.5G/DL HI: 6.0 G/DL
NORMAL RANGE:	LO: 3.5 G/DL HI: 5.3 G/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR CAL-1:	4.0 G/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.600 HI: 0.800
BLANK RANGE:	LO: 0.060 HI: 1.100

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: PHOSPHORUS
Catalog # : P7516

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

R1 REAGENT 1:	200 UL
I1 INCUBATION:	5 S
M1:	4 S
S SAMPLE:	8 UL
DILUENT:	30 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	90 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 12.0
NORMAL RANGE:	LO: 2.5 HI: 4.8 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH TEST
CALIBRATOR CAL-1:	5.0 MG/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: NO HI: NO
BLANK RANGE:	LO: NO HI: NO

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test: TOTAL PROTEIN

Catalog # : T7528

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	R1-I1-S-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	540 NM
TEMPERATURE:	37C
DECIMAL POSITION:	2
UNIT:	G/DL

Analysis:

R1 REAGENT 1:	200
I1 INCUBATION:	5 S
M1:	4 S
S SAMPLE:	4 UL
DILUENT:	30 UL

AO READINGS

FIRST:	0.5 S
NUMBER:	2
INTERVAL:	115 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 1 HI: 15.0 G/DL
NORMAL RANGE:	LO: 6.2 HI: 8.5 G/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2
REAC. LIMIT:	
POINT:	

Calibration:

CALIB. INTERVAL:	EACH RUN
CALIBRATOR CAL-1:	8.0 G/DL
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.050 HI: 0.200
BLANK RANGE:	LO: -0.001 HI: 0.010

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara

Test:

Catalog # :

General:

MEASUREMENT MODE:
REACTION MODE:
CALIBRATION MODE:
REAGENT BLANK:
WAVELENGTH:
TEMPERATURE:
DECIMAL POSITION:
UNIT:

Analysis:

R1 REAGENT 1:
I1 INCUBATION:
M1:
S SAMPLE:
DILUENT:

AO READINGS
FIRST:
NUMBER:
INTERVAL:

Calculation:

CONVERS. FACTOR:
OFFS:
REAC. DIRECTION:
CHECK:
SAMPLE LIMIT:
TEST RANGE: LO: HI:
NORMAL RANGE: LO: HI:
CALC. STEPS:
CALC. STEP A:
READINGS: FIRST: LAST:
REAC. LIMIT:
POINT:

Calibration:

CALIB. INTERVAL:
CALIBRATOR CAL-1:
POS:
REPLICATE:
REAGENT RANGE: LO: HI:
BLANK RANGE: LO: HI:

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: DIRECT BILIRUBIN
Catalog # : B7538

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	546 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	180 S
I1 INCUBATION:	185 S
P SAMPLE:	20 UL
DILUENT:	30 UL
SR 1 START R1:	20 UL
REAGENT:	100 UL
DILUENT:	30 UL
AO READINGS	
FIRST:	1.0 S
NUMBER:	2
INTERVAL:	300 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 0.0 HI: 0.5 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	269.0
POS:	
REPLICATE:	
REAGENT RANGE:	LO: -0.0150 HI: 0.0050
BLANK RANGE:	LO: -0.0050 HI: 0.015

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: TOTAL BILIRUBIN
Catalog # : B7576

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	FACTOR
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	546 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	10 S
I1 INCUBATION:	15 S
P SAMPLE:	20 UL
DILUENT:	30 UL
SR 1 START R1:	20 UL
REAGENT:	100 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	1.0 S
NUMBER:	2
INTERVAL:	180 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 0.2 HI: 1.2 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	
FACTOR:	269.800
POS:	
REPLICATE:	
REAGENT RANGE:	LO: -0.0150 HI: 0.0050
BLANK RANGE:	LO: -0.0050 HI: 0.015

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CALCIUM
Catalog # : C7503, C7508

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	573 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MG/DL

Analysis:

M1:	35 S
I1 INCUBATION:	40 S
P SAMPLE:	10 UL
DILUENT:	50 UL
SR 1 START R1:	50 UL
REAGENT:	150 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	0.5 S
NUMBER:	2
INTERVAL:	40 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 0.0 HI: 20.0 MG/DL
NORMAL RANGE:	LO: 8.5 HI: 10.4 MG/DL
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 2

Calibration

CALIB. INTERVAL:	ON REQUEST
CALIBRATOR:	10.0 MG/DL
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: -0.100 HI: 0.500
BLANK RANGE:	LO: 0.100 HI: 0.400

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: CHLORIDE
Catalog # : C7501

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-AO
CALIBRATION MODE:	CALIB
REAGENT BLANK:	REAG/DIL
WAVELENGTH:	480 NM
TEMPERATURE:	37C
DECIMAL POSITION:	1
UNIT:	MEQ/L

Analysis:

M1:	120 S
I1 INCUBATION:	125 S
P SAMPLE:	5 UL
DILUENT:	25 UL
SR 1 START R1:	50 UL
REAGENT:	150 UL
DILUENT:	20 UL
AO READINGS	
FIRST:	60.0 S
NUMBER:	1
INTERVAL:	5 S

Calculation:

CONVERS. FACTOR:	1.000
OFFS:	0.000
REAC. DIRECTION:	INCREASE
CHECK:	ON
SAMPLE LIMIT:	NO
TEST RANGE:	LO: 80 HI: 120 MEQ/L
NORMAL RANGE:	LO: 98 HI: 106 MEQ/L
CALC. STEPS:	
CALC. STEP A:	
READINGS:	FIRST: M1 LAST: 1

Calibration

CALIB. INTERVAL:	EACH TEST
CALIBRATOR:	100 MEQ/L
FACTOR:	
POS:	3
REPLICATE:	SINGLE
REAGENT RANGE:	LO: 0.0000 HI: 0.3000
BLANK RANGE:	LO: -0.010 HI: 0.002

Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Fara
Test: Lp(a)
Catalog # : L7597

General:

MEASUREMENT MODE:	ABSORB
REACTION MODE:	P-I1-SR1-A0
CALIBRATION MODE:	STD Non Linear
REAGENT BLANK:	RGT/DIL
WAVELENGTH:	340 NM
TEMPERATURE:	37C
DECIMAL POSITION:	0
UNIT:	MG/DL

Analysis:

P SAMPLE:	11UL
DILUENT:	30 UL
REAGENT:	210 UL
I1 Incubation	240 S
M1	220 S
SR1 Start Reagent 1	35 S
Diluent	20 S
AO READINGS	
FIRST:	5 S
NUMBER:	6 S
INTERVAL:	240 S
M1:	

Calculation:

CONVERS. FACTOR:	1.000	
OFFS:	0.000	
REAC. DIRECTION:	INCREASE	
CHECK:	OFF	
SAMPLE LIMIT:	NO	
TEST RANGE:	LO: NO	HI:NO
NORMAL RANGE:	LO: 0	HI:30
CALC. STEPS:	1	
CALC. STEP A:	Endpoint	
READINGS:	FIRST: M1	LAST: 6
REAC. LIMIT:		
POINT:		

Calibration:

CALIB. INTERVAL: On Request	Replicate : 2 Duplicate
Deviation %: 10	Calculation Model : Linear Interpolation

Calibration Correction : NO

Antigen Excess : NO

