

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917  
 Test: CRP  
 Catalog # :C7564

Rev. 12-02

## ANALYZE

TEST NAME	[CRP]			
ASSAY/ POINT	[2-POINT END]	[10]	[19]	[34][0][0]
WAVE (SUB/MAIN)	[800]	[570]		
S. VOL. (NORMAL)	[15.0]	[0.0]	[0.0]	
S. VOL. (DECREASE)	[15.0]	[0.0]	[0.0]	
S. VOL. (INCREASE)	[15.0]	[0.0]	[0.0]	
DILUENT	[DW]	[0]		
REAGENT VOL (R1)	[150]	[0]	[ ]	[0]
REAGENT VOL (R2)	[0]	[0]	[ ]	[0]
REAGENT VOL (R3)	[100]	[0]	[ ]	[0]
REAGENT VOL (R4)	[0]	[0]	[ ]	[0]
ABS. LIMIT	[32000]	[0]	[INCREASE]	TWIN: [ ]
PROZONE LIMIT:	[-32000]	[0]	[LOWER]	
CELL DETERGENT	[DETERGENT]			

## CALIBRATION

CALIB TYPE	[SPLINE]	[ ]	
POINT	[6]	SPAN	POINTE [6]
WEIGHT	[0]		
AUTO TIME OUTN	Auto change		
BLANK	[0]	CHANGE LOT:	[ ]
SPAN	[0]	CHANGE BOTTLE:	[ ]
2 POINT	[0]		
FULL	[0]		
SD LIMIT	[999.9]	DUPLICATE:	[32000]
SENSITIVITY	[0]	S1ABS RANGE:	[-32000] [32000]

## RANGE

TEST CODE	[ ]	UNITS	:[MG/DL]	DATA MODE	:[ON BOARD]
CONTROL INTERVAL	[0]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[-999999]		[999999]		

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>						
CONC	[0.00]	[0.05]	[0.15]	[0.50]	[0.75]	[1.0]
POSITION	[*]	[*]	[*]	[*]	[*]	[*]
SAMPLE :	[15.0]	[15.0]	[15.0]	[15.0]	[15.0]	[15.0]
<u>PRE-DILUENT</u>						
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
CALIB CODE:	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]

Use saline solution at blank (Std1)

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917  
 Test: LDL-CHOLESTEROL  
 Catalog # :L7574

Rev. 12-02

## ANALYZE

TEST NAME	[LDL]				
ASSAY/ POINT	[2-POINT END]	[10]	[16]	[34]	[0] [0]
WAVE (SUB/MAIN)	[700]	[546]			
S. VOL. (NORMAL)	[3.0]	[0.0]	[0.0]		
S. VOL. (DECREASE)	[2.0]	[0.0]	[0.0]		
S. VOL. (INCREASE)	[3.0]	[0.0]	[0.0]		
DILUENT	[00951]	[99]			
REAGENT VOL (R1)	[200]	[0]	[00410]	[28]	
REAGENT VOL (R2)	[0]	[0]	[00410]	[0]	
REAGENT VOL (R3)	[67]	[0]	[00410]	[28]	
REAGENT VOL (R4)	[0]	[0]	[00410]	[ ]	
ABS. LIMIT	[ ]	[0]	[INCREASE]		TWIN: [ ]
PROZONE LIMIT:	[32000]	[0]	[UPPER]		
CELL DETERGENT	[DETERGENT]				

## CALIBRATION

CALIB TYPE	[LINEAR]	[ ]		
POINT	[2]	SPAN	POINTE	[2]
WEIGHT	[0]			
AUTO TIME OUTN	Auto change			
BLANK	[0]	CHANGE LOT:	[ ]	
SPAN	[0]	CHANGE BOTTLE:	[ ]	
2 POINT	[0]			
FULL	[0]			
SD LIMIT	[0.1]	DUPLICATE:	[5%]	
SENSITIVITY	[10]	S1ABS RANGE:	[ ] [ ]	

## RANGE

TEST CODE	[449]	UNITS	:[MG/DL]	DATA MODE	:[ON BOARD]
CONTROL INTERVAL	[1000]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[-999999]		[999999]		

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>						
CONC	[501]	[674]	[0]	[0]	[0]	[0]
POSITION	[*]	[*]	[*]	[*]	[*]	[*]
SAMPLE :	[3.0]	[3.0]	[0.0]	[0.0]	[0.0]	[0.0]
<u>PRE-DILUENT</u>						
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0]	[0]	[0]	[0]	[0]	[0]
CALIB CODE:	[ ]	[ ]	[ ]	[ ]	[ ]	[ ]

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917

Test: Iron

Catalog # :I7504

Rev: 1-03 non validated

application

## ANALYZE

TEST NAME	[Iron]				
ASSAY/ POINT	[2-POINT END]	[10]	[16]	[20]	[00][0]
WAVE (SUB/MAIN)	[700]	[570]			
S. VOL. (NORMAL)	[15.0]	[0.0]	[0.0]		
S. VOL. (DECREASE)	[25.0]	[0.0]	[0.0]		
S. VOL. (INCREASE)	[15.0]	[0.0]	[0.0]		
DILUENT	[ ]	[99]			
REAGENT VOL (R1)	[180]	[0]	[ ]	[0]	
REAGENT VOL (R2)	[0]	[0]	[ ]	[0]	
REAGENT VOL (R3)	[36]	[0]	[ ]	[0]	
REAGENT VOL (R4)	[0]	[0]	[ ]	[ ]	
ABS. LIMIT	[32000]	[0]	[INCREASE]		TWIN: [ ]
PROZONE LIMIT:	[-32000]	[0]	[LOWER]		
CELL DETERGENT	[DETERGENT]				

## CALIBRATION

CALIB TYPE	[LINEAR]	[ ]		
POINT	[2]	SPAN POINTE	[2]	
WEIGHT	[0]			
AUTO TIME OUTN	Auto change			
BLANK	[24]	CHANGE LOT:	[2pointA]	
SPAN	[0]	CHANGE BOTTLE:	[ ]	
2 POINT	[blankA]			
FULL	[0]			
SD LIMIT	[0.1]	DUPLICATE:	[32000]	
SENSITIVITY	[0]	S1ABS RANGE:	[0]	[4000]

## RANGE

TEST CODE	[ ]	UNITS	:[ug/dl]	DATA MODE	:[ON BOARD]
CONTROL INTERVAL	[0]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[-999999]		[999999]		

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>					
CONC	[0.00]	[*]	[0]	[0]	[0]
POSITION	[*]	[*]	[ ]	[ ]	[ ]
SAMPLE :	[15.0]	[15.0]	[15.0]	[15.0]	[15.0]
<u>PRE-DILUENT</u>					
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
CALIB CODE:	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]

Use saline solution at blank (Std1)

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.



# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917  
 Test: RF  
 Catalog # :R7568

rev: 2-03 non validated

application

## ANALYZE

TEST NAME	[RF ]				
ASSAY/ POINT	[2-POINT END ]	[10 ]	[16 ]	[34 ]	[0 ] [0 ]
WAVE (SUB/MAIN)	[700 ]	[340 ]			
S. VOL. (NORMAL)	[12.0 ]	[0.0 ]	[0.0 ]		
S. VOL. (DECREASE)	[0.0 ]	[0.0 ]	[0.0 ]		
S. VOL. (INCREASE)	[20.0 ]	[0.0 ]	[0.0 ]		
DILUENT	[DW ]	[0 ]			
REAGENT VOL (R1)	[200 ]	[0 ]	[0 ]	[0 ]	
REAGENT VOL (R2)	[0 ]	[0 ]	[0 ]	[0 ]	
REAGENT VOL (R3)	[60 ]	[0 ]	[0 ]	[0 ]	
REAGENT VOL (R4)	[0 ]	[0 ]	[0 ]	[0 ]	
ABS. LIMIT	[ ]	[ ]	[ ]		TWIN: [ ]
PROZONE LIMIT:	[0 ]	[0 ]	[LOWER ]		
CELL DETERGENT	[ ]				

## CALIBRATION

CALIB TYPE	[SPLINE ]	[ ]		
POINT	[6 ]	SPAN POINTE	[5 ]	
WEIGHT	[0 ]			
AUTO TIME OUTN	Auto change			
BLANK	[0 ]	CHANGE LOT:	[ ]	
SPAN	[0 ]	CHANGE BOTTLE:	[ ]	
2 POINT	[ ]			
FULL	[0 ]			
SD LIMIT	[999.9 ]	DUPLICATE:	[500 ]	
SENSITIVITY	[0 ]	S1ABS RANGE:	[-32000 ] [32000 ]	

## RANGE

TEST CODE	[ ]	UNITS	:[U/ml ]	DATA MODE	:[ ]
CONTROL INTERVAL	[ ]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0] ]		
PANIC VALUE	[ ]	[ ]			

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>					
CONC	[0.0 ]	[* ]	[* ]	[* ]	[ ] [ ]
POSITION	[1 ]	[* ]	[* ]	[* ]	[ ] [ ]
SAMPLE :	[12.0 ]	[12.0 ]	[12.0 ]	[12.0 ]	[ ] [ ]
<u>PRE-DILUENT</u>					
VOLUME	[0.0 ]	[0.0 ]	[0.0 ]	[0.0 ]	[0.0 ] [0.0 ]
DILUENT	[0 ]	[0 ]	[0 ]	[0 ]	[0 ] [0 ]
CALIB CODE:	[0 ]	[0 ]	[0 ]	[0 ]	[0 ] [0 ]

\* User Defined

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917  
 Test: ASO  
 Catalog # :A7566

rev: 2-03 non validated

application

## ANALYZE

TEST NAME	[ASO]				
ASSAY/ POINT	[2-POINT END]	[10]	[16]	[34]	[0] [0]
WAVE (SUB/MAIN)	[700]	[340]			
S. VOL. (NORMAL)	[15.0]	[0.0]	[0.0]		
S. VOL. (DECREASE)	[10.0]	[0.0]	[0.0]		
S. VOL. (INCREASE)	[30.0]	[0.0]	[0.0]		
DILUENT	[DW]	[0]			
REAGENT VOL (R1)	[300]	[0]	[0]	[0]	
REAGENT VOL (R2)	[0]	[0]	[0]	[0]	
REAGENT VOL (R3)	[100]	[0]	[0]	[0]	
REAGENT VOL (R4)	[0]	[0]	[0]	[0]	
ABS. LIMIT	[ ]	[ ]	[ ]	TWIN:	[ ]
PROZONE LIMIT:	[-32000]	[34]	[LOWER]		
CELL DETERGENT	[ ]				

## CALIBRATION

CALIB TYPE	[SPLINE]	[ ]		
POINT	[6]	SPAN POINTE	[6]	
WEIGHT	[0]			
AUTO TIME OUTN	Auto change			
BLANK	[0]	CHANGE LOT:	[ ]	
SPAN	[0]	CHANGE BOTTLE:	[ ]	
2 POINT	[ ]			
FULL	[0]			
SD LIMIT	[999]	DUPLICATE:	[10000]	
SENSITIVITY	[0]	S1ABS RANGE:	[-32000] [32000]	

## RANGE

TEST CODE	[ ]	UNITS	:[U/ml]	DATA MODE	:[ ]
CONTROL INTERVAL	[ ]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[ ]	[ ]			

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>					
CONC	[0.0]	[*]	[*]	[*]	[ ] [ ]
POSITION	[1]	[*]	[*]	[*]	[ ] [ ]
SAMPLE :	[15.0]	[15.0]	[15.0]	[15.0]	[ ] [ ]
<u>PRE-DILUENT</u>					
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0] [0.0]
DILUENT	[0]	[0]	[0]	[0]	[0] [0]
CALIB CODE:	[0]	[0]	[0]	[0]	[0] [0]

\* User Defined

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917

Test: FRU

Catalog # :F7546

rev: 3-03

## ANALYZE

TEST NAME	[Fructosamine ]				
ASSAY/ POINT	[2-POINT END ]	[10]	[31]	[34]	[0] [0]
WAVE (SUB/MAIN)	[700]	[546]			
S. VOL. (NORMAL)	[20.0]	[0.0]	[0.0]		
S. VOL. (DECREASE)	[20.0]	[0.0]	[0.0]		
S. VOL. (INCREASE)	[20.0]	[0.0]	[0.0]		
DILUENT	[DW]	[0]			
REAGENT VOL (R1)	[250]	[0]	[0]	[0]	
REAGENT VOL (R2)	[0]	[0]	[0]	[0]	
REAGENT VOL (R3)	[0]	[0]	[0]	[0]	
REAGENT VOL (R4)	[0]	[0]	[0]	[0]	
ABS. LIMIT	[32000]	[0]	[INCREASE]		TWIN: []
PROZONE LIMIT:	[32000]	[0]	[LOWER]		
CELL DETERGENT	[DETERGENT]				

## CALIBRATION

CALIB TYPE	[LINEAR]	[ ]		
POINT	[2]	SPAN	POINTE	[2]
WEIGHT	[0]			
AUTO TIME OUTN	Auto change			
BLANK	[0]	CHANGE LOT:	[ ]	
SPAN	[0]	CHANGE BOTTLE:	[ ]	
2 POINT	[4]			
FULL	[0]			
SD LIMIT	[0.1]	DUPLICATE:	[500]	
SENSITIVITY	[0]	S1ABS RANGE:	[-32000] [32000]	

## RANGE

TEST CODE	[ ]	UNITS	:[mmol/L]	DATA MODE	:[ON BOARD]
CONTROL INTERVAL	[0]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[-99999]	[99999]			

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>						
CONC	[0.0]	[2.9]	[ ]	[ ]	[ ]	[ ]
POSITION	[*]	[*]	[ ]	[ ]	[ ]	[ ]
SAMPLE :	[20.0]	[20.0]	[ ]	[ ]	[ ]	[ ]
<u>PRE-DILUENT</u>						
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0]	[0]	[0]	[0]	[0]	[0]
CALIB CODE:	[0]	[0]	[0]	[0]	[0]	[0]

\* User Defined

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917

Test: BHY

Catalog # :H7587

## ANALYZE

TEST NAME	[Beta Hydroxy ]			
ASSAY/ POINT	[2-POINT END ]	[10]	[16]	[34] [0] [0]
WAVE (SUB/MAIN)	[-]	[505]		
S. VOL. (NORMAL)	[8.0]	[0.0]	[0.0]	
S. VOL. (DECREASE)	[0.0]	[0.0]	[0.0]	
S. VOL. (INCREASE)	[0.0]	[0.0]	[0.0]	
DILUENT	[DW]	[0]		
REAGENT VOL (R1)	[300]	[0]	[0]	[0]
REAGENT VOL (R2)	[0]	[0]	[0]	[0]
REAGENT VOL (R3)	[50]	[0]	[0]	[0]
REAGENT VOL (R4)	[0]	[0]	[0]	[0]
ABS. LIMIT	[-32000]	[0]	[INCREASE]	TWIN: []
PROZONE LIMIT:	[0]	[0]	[LOWER]	
CELL DETERGENT	[Detergent ]			

## CALIBRATION

CALIB TYPE	[LINEAR]	[ ]
POINT	[2]	SPAN POINTE [2]
WEIGHT	[0]	
AUTO TIME OUTN	Auto change	
BLANK	[0]	CHANGE LOT: [0]
SPAN	[0]	CHANGE BOTTLE: [ ]
2 POINT	[0]	
FULL	[0]	
SD LIMIT	[0.1]	DUPLICATE: [100]
SENSITIVITY	[0]	S1ABS RANGE: [-32000] [32000]

## RANGE

TEST CODE	[*]	UNITS :[mM]	DATA MODE :[ON BOARD]
CONTROL INTERVAL	[1000]	INST. FACTOR: [(Y=aX=b) : a=[1.0] b = [0.0]]	
PANIC VALUE	[-99999]	[99999]	

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>						
CONC	[0.0]	[1.0]	[0]	[0]	[0]	[0]
POSITION	[*]	[*]	[ ]	[ ]	[ ]	[ ]
SAMPLE :	[8.0]	[8.0]	[ ]	[ ]	[ ]	[ ]
<u>PRE-DILUENT</u>						
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0.0]	[0.]	[0.0]	[0.0]	[0.0]	[0.0]
CALIB CODE:	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]

\*Non Validated application



# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 917

Test: HbA1c

Catalog # :H7541

rev: 8-03

## ANALYZE

TEST NAME	[HbA1c ]			
ASSAY/ POINT	[2-POINT END ]	[10]	[19]	[34] [0] [0]
WAVE (SUB/MAIN)	[none ]	[660]		
S. VOL. (NORMAL)	[4.0]	[0.0]	[0.0]	
S. VOL. (DECREASE)	[4.0]	[0.0]	[0.0]	
S. VOL. (INCREASE)	[4.0]	[0.0]	[0.0]	
DILUENT	[DW]	[0]		
REAGENT VOL (R1)	[150]	[0]	[0]	[0]
REAGENT VOL (R2)	[0]	[0]	[0]	[0]
REAGENT VOL (R3)	[50]	[0]	[0]	[0]
REAGENT VOL (R4)	[0]	[0]	[0]	[0]
ABS. LIMIT	[25000]	[0]	[INCREASE]	TWIN: []
PROZONE LIMIT:	[-32000]	[0]	[LOWER]	
CELL DETERGENT	[Detergent]			

## CALIBRATION

CALIB TYPE	[SPLINE]	[ ]		
POINT	[5]	SPAN	POINTE	[5]
WEIGHT	[0]			
AUTO TIME OUTN	Auto change			
BLANK	[0]	CHANGE LOT:	[0]	
SPAN	[0]	CHANGE BOTTLE:	[ ]	
2 POINT	[0]			
FULL	[0]			
SD LIMIT	[999]	DUPLICATE:	[32000]	
SENSITIVITY	[0]	S1ABS RANGE:	[-32000] [32000]	

## RANGE

TEST CODE	[ ]	UNITS	:[%]	DATA MODE	:[ON BOARD]
CONTROL INTERVAL	[ ]	INST. FACTOR:	[(Y=aX=b) : a=[1.0] b = [0.0]]		
PANIC VALUE	[-99999]	[99999]			

## STANDARD CONCENTRATION

\*USER DEFINED

<u>STANDARD SOLUTION</u>						
CONC	[0.0]	[*]	[*]	[*]	[*]	[ ]
POSITION	[*]	[*]	[*]	[*]	[*]	[ ]
SAMPLE :	[4.0]	[4.0]	[4.0]	[4.0]	[4.0]	[ ]
<u>PRE-DILUENT</u>						
VOLUME	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
DILUENT	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]
CALIB CODE:	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]	[0.0]

\* User Defined

It is recommended that two levels of control material be assayed daily.

