

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E

Test: HbA1c

Catalog #: H7541

TEST PARAMETERS

Name	:HbA1c	Ref. Male low	:0.00%
Abbr. Name	:HbA1c	Ref. Male high	:6.00%
Mode	:Twopoint	Ref. Female low	:0.00%
Wavelength	:620 nm	Ref. Female high	:6.00%
Units	:%	Ref. Ped. Low	:0.00%
Decimals	:1	Ref. Ped. High	:6.00%
Low Conc.	:2.00%	Control 1	:User Defined
High Conc.	:16.00%	Control 2	:User Defined
Calibrator Name	:HbA1c Cal	Control 3	:User Defined
Repeat	:1	Correlat. Factor	:1.000
Number	:4	Correlat offset	:0.000
Interval	:1 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:180 ul
Rerun volume	:180 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
R2 bottle	:5 ml
Normal volume	:60 ul
Rerun volume	:60 ul
Predilution	:No
Slope Blank	:No
Point one, two	:6, 236 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:2.000
Reagent Blank	:No
Max. inaccuracy	:1.00%
Prozone Check	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Mono Mode not applicable

It is recommended that two levels of control material be assayed daily.

Rev 8-03 non validated application

Pointe Scientific Test Parameters - Selectra E

Total Bilirubin (TBILI): B7576

TEST PARAMETERS

Name	:Total Bilirubin	Prozone Check	:No
Abbr. Name	:TBILI	Ref. Male low	:0.2 mg/dl
Mode	:Endpoint	Ref. Male high	:1.2 mg/dl
Wavelength	:546 nm	Ref. Female low	:0.2 mg/dl
Units	:mg/dl	Ref. Female high	:1.2 mg/dl
Decimals	:1	Ref. Ped. Low	:0.2 mg/dl
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:1.2 mg/dl
High Conc.	:20.0 mg/dl	Control 1	:User Defined
Calibrator Name	:Bili Cal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.800
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.800
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Calcium (CA): C7529

TEST PARAMETERS

Name	:Calcium	Prozone Check	:No
Abbr. Name	:CA	Ref. Male low	:8.5 mg/dl
Mode	:Endpoint	Ref. Male high	:10.4 mg/dl
Wavelength	:620 nm	Ref. Female low	:8.5 mg/dl
Units	:mg/dl	Ref. Female high	:10.4 mg/dl
Decimals	:1	Ref. Ped. Low	:8.5 mg/dl
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:10.4 mg/dl
High Conc.	:20.0 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:4.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:2.500
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:2.500
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

autoHDL (HDLC): H7545

TEST PARAMETERS

Name	:autoHDL	Prozone Check	:No
Abbr. Name	:HDLC	Ref. Male low	:35 mg/dl
Mode	:Endpoint	Ref. Male high	:65 mg/dl
Wavelength	:578 nm	Ref. Female low	:35 mg/dl
Units	:mg/dl	Ref. Female high	:65 mg/dl
Decimals	:0	Ref. Ped. Low	:35 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:65 mg/dl
High Conc.	:150 mg/dl	Control 1	:User Defined
Calibrator Name	:autoHDL Cal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:295 ul
Rerun volume	:295 ul
Sample	
Normal volume	:4.0 ul
Rerun volume	:4.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:100 ul
Rerun volume	:100 ul
Predilution	:No
Incubation time	:6.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:
R1 bottle (ml)	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Incubation time	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
Reagent Blank	:
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

autoLDL (LDLC): L7574

TEST PARAMETERS

Name	:autoLDL	Prozone Check	:No
Abbr. Name	:LDLC	Ref. Male low	:0 mg/dl
Mode	:Endpoint	Ref. Male high	:130 mg/dl
Wavelength	:546 nm	Ref. Female low	:0 mg/dl
Units	:mg/dl	Ref. Female high	:130 mg/dl
Decimals	:0	Ref. Ped. Low	:0 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:130 mg/dl
High Conc.	:700 mg/dl	Control 1	:User Defined
Calibrator Name	:autoLDL Cal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:295 ul
Rerun volume	:295 ul
Sample	
Normal volume	:4.0 ul
Rerun volume	:4.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:100 ul
Rerun volume	:100 ul
Predilution	:No
Incubation time	:6.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:
R1 bottle (ml)	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Incubation time	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
Reagent Blank	:
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

Cholesterol (CHOL): C7510

TEST PARAMETERS

Name	:Cholesterol	Prozone Check	:No
Abbr. Name	:CHOL	Ref. Male low	:0 mg/dl
Mode	:Endpoint	Ref. Male high	:200 mg/dl
Wavelength	:505-620	Ref. Female low	:0 mg/dl
Units	:mg/dl	Ref. Female high	:200 mg/dl
Decimals	:0	Ref. Ped. Low	:0 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:200 mg/dl
High Conc.	:700 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.200
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.200
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Glucose (HEX) (GLU): G7517

TEST PARAMETERS

Name	:Glucose (HEX)	Prozone Check	:No
Abbr. Name	:GLU	Ref. Male low	:65 mg/dl
Mode	:Endpoint	Ref. Male high	:110 mg/dl
Wavelength	:340-620	Ref. Female low	:65 mg/dl
Units	:mg/dl	Ref. Female high	:110 mg/dl
Decimals	:0	Ref. Ped. Low	:65 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:110 mg/dl
High Conc.	:600 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.500
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:11.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.500
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Iron (FE): 17505

TEST PARAMETERS

Name	:Iron	Prozone Check	:No
Abbr. Name	:FE	Ref. Male low	:60 ug/dl
Mode	:Endpoint	Ref. Male high	:150 ug/dl
Wavelength	:546 nm	Ref. Female low	:60 ug/dl
Units	:ug/dl	Ref. Female high	:150 ug/dl
Decimals	:0	Ref. Ped. Low	:60 ug/dl
Low Conc.	:0 ug/dl	Ref. Ped. High	:150 ug/dl
High Conc.	:900 ug/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:Yes
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:25.0 ul
Rerun volume	:25.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:25.0 ul
Rerun volume	:25.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Phosphorus (PHOS): P7516

TEST PARAMETERS

Name	:Phosphorus	Prozone Check	:No
Abbr. Name	:PHOS	Ref. Male low	:2.5 mg/dl
Mode	:Endpoint	Ref. Male high	:4.8 mg/dl
Wavelength	:340-620 nm	Ref. Female low	:2.5 mg/dl
Units	:mg/dl	Ref. Female high	:4.8 mg/dl
Decimals	:1	Ref. Ped. Low	:2.5 mg/dl
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:4.8 mg/dl
High Conc.	:20.0 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.800
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.800
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Total Protein (TP): T7528

TEST PARAMETERS

Name	:Total Protein	Prozone Check	:No
Abbr. Name	:TP	Ref. Male low	:6.2 g/dl
Mode	:Endpoint	Ref. Male high	:8.5 g/dl
Wavelength	:546 nm	Ref. Female low	:6.2 g/dl
Units	:g/dl	Ref. Female high	:8.5 g/dl
Decimals	:1	Ref. Ped. Low	:6.2 g/dl
Low Conc.	:0.0 g/dl	Ref. Ped. High	:8.5 g/dl
High Conc.	:15.0 g/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.300
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Triglyceride (GPO) (TRIG): T7532

TEST PARAMETERS

Name	:Triglyceride (GPO)	Prozone Check	:No
Abbr. Name	:TRIG	Ref. Male low	:36 mg/dl
Mode	:Endpoint	Ref. Male high	:165 mg/dl
Wavelength	:505 nm	Ref. Female low	:36 mg/dl
Units	:mg/dl	Ref. Female high	:165 mg/dl
Decimals	:0	Ref. Ped. Low	:36 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:165 mg/dl
High Conc.	:1000 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:1.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:1.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Uric Acid (UA): U7581

TEST PARAMETERS

Name	:Uric Acid	Prozone Check	:No
Abbr. Name	:UA	Ref. Male low	:2.5 mg/dl
Mode	:Endpoint	Ref. Male high	:7.7 mg/dl
Wavelength	:505-620 nm	Ref. Female low	:2.5 mg/dl
Units	:mg/dl	Ref. Female high	:7.7 mg/dl
Decimals	:1	Ref. Ped. Low	:2.5 mg/dl
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:7.7 mg/dl
High Conc.	:25.0 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:8.0 ul
Rerun volume	:8.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.600
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:8.0 ul
Rerun volume	:8.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.600
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Magnesium (MG): M7527

TEST PARAMETERS

Name	:Magnesium	Prozone Check	:No
Abbr. Name	:MG	Ref. Male low	:1.3 mEq/l
Mode	:Endpoint	Ref. Male high	:2.5 mEq/l
Wavelength	:546 nm	Ref. Female low	:1.3 mEq/l
Units	:mEq/l	Ref. Female high	:2.5 mEq/l
Decimals	:1	Ref. Ped. Low	:1.3 mEq/l
Low Conc.	:0.0 mEq/l	Ref. Ped. High	:2.5 mEq/l
High Conc.	:6.0 mEq/l	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:4.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Protein, Urine (UPRO): P7582

TEST PARAMETERS

Name	:Protein, Urine	Prozone Check	:No
Abbr. Name	:UPRO	Ref. Male low	:28 mg/24h
Mode	:Endpoint	Ref. Male high	:141 mg/24h
Wavelength	:620 nm	Ref. Female low	:28 mg/24h
Units	:mg/dl	Ref. Female high	:141 mg/24h
Decimals	:1	Ref. Ped. Low	:28 mg/24h
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:141 mg/24h
High Conc.	:250.0 mg/dl	Control 1	:User Defined
Calibrator Name	:UPRO STD	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:6.0 ul
Rerun volume	:6.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:4.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:6.0 ul
Rerun volume	:6.0 ul
Incubation time	:11.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Ammonia (AMM): A7553

TEST PARAMETERS

Name	:Ammonia	Prozone Check	:No
Abbr. Name	:AMM	Ref. Male low	:11 umol/l
Mode	:Endpoint	Ref. Male high	:35 umol/l
Wavelength	:340 nm	Ref. Female low	:11 umol/l
Units	:umol/l	Ref. Female high	:35 umol/l
Decimals	:0	Ref. Ped. Low	:11 umol/l
Low Conc.	:0 umol/l	Ref. Ped. High	:35 umol/l
High Conc.	:600 umol/l	Control 1	:User Defined
Calibrator Name	:AMM STD	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:2	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:240 ul
Rerun volume	:240 ul
Sample	
Normal volume	:30.0 ul
Rerun volume	:30.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:13 ul
Rerun volume	:13 ul
Predilution	:No
Incubation time	:6.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.800
R. Abs. H. Limit	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:
R1 bottle (ml)	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Incubation time	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
Reagent Blank	:
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

Ethanol (ETOH): A7504

TEST PARAMETERS

Name	:Ethanol	Prozone Check	:No
Abbr. Name	:ETOH	Ref. Male low	:0 mg/dl
Mode	:Endpoint	Ref. Male high	:100 mg/dl
Wavelength	:340 nm	Ref. Female low	:0 mg/dl
Units	:mg/dl	Ref. Female high	:100 mg/dl
Decimals	:0	Ref. Ped. Low	:0 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:100 mg/dl
High Conc.	:400 mg/dl	Control 1	:User Defined
Calibrator Name	:ETOH STD	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:395 ul
Rerun volume	:395 ul
Sample	
Normal volume	:2.0 ul
Rerun volume	:2.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:4.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:395 ul
Rerun volume	:395 ul
Sample	
Normal volume	:2.0 ul
Rerun volume	:2.0 ul
Incubation time	:11.5 min
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Lactic Acid (LAC): L7596

TEST PARAMETERS

Name	:Lactic Acid	Prozone Check	:No
Abbr. Name	:LAC	Ref. Male low	:0.5 mmol/l
Mode	:Endpoint	Ref. Male high	:2.2 mmol/l
Wavelength	:546 nm	Ref. Female low	:0.5 mmol/l
Units	:mmol/l	Ref. Female high	:2.2 mmol/l
Decimals	:1	Ref. Ped. Low	:0.5 mmol/l
Low Conc.	:0.0 mmol/l	Ref. Ped. High	:2.2 mmol/l
High Conc.	:20.0 mmol/l	Control 1	:User Defined
Calibrator Name	:LAC STD	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:30 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:225 ul
Rerun volume	:225 ul
Sample	
Normal volume	:4.0 ul
Rerun volume	:4.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:150 ul
Rerun volume	:150 ul
Predilution	:No
Incubation time	:6.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:375 ul
Rerun volume	:375 ul
Sample	
Normal volume	:4.0 ul
Rerun volume	:4.0 ul
Incubation time	:11.5 min
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Carbon Dioxide (CO2): C7504

TEST PARAMETERS

Name	:Carbon Dioxide	Prozone Check	:No
Abbr. Name	:CO2	Ref. Male low	:23 mmol/l
Mode	:Twopoint	Ref. Male high	:34 mmol/l
Wavelength	:340 nm	Ref. Female low	:23 mmol/l
Units	:mmol/l	Ref. Female high	:34 mmol/l
Decimals	:0	Ref. Ped. Low	:23 mmol/l
Low Conc.	:0.0 mmol/l	Ref. Ped. High	:34 mmol/l
High Conc.	:40 mmol/l	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Incubation time	:24, 77 sec
Low absorbance	:0.300
High absorbance	:3.000
R. Abs. L. Limit	:0.700
R. Abs. H. Limit	:3.000
Reagent Blank	:Yes
Cal. High Limit	:0.000 Abs
Cal. Low Limit	:0.000 Abs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Incubation time	:12, 70 sec
Low absorbance	:0.300
High absorbance	:3.000
R. Abs. L. Limit	:0.700
R. Abs. H. Limit	:3.000
Reagent Blank	:No
Cal. Low Limit	:0.000 Abs
Cal. High Limit	:0.000 Abs
Factor	:1.00

NOTES:

REV:3/00



Analyzer: Selectra E
Test: Acid Phos
Catalog # A7503

TEST PARAMETERS

Name	:Acid Phosphatase	Ref. Ped. Low	:0 U/I
Abbr. Name	:ACP	Ref. Ped. High	:11 U/I
Mode	:Kinetic	Control 1	:User Defined
Wavelength	:405 nm	Control 2	:User Defined
Units	:U/I	Control 3	:User Defined
Decimals	:0	Correlat. Factor	:1.000
Low Conc.	:0 U/I	Correlat offset	:0.000
High Conc.	:35 U/I		
Calibrator Name	:none		
Prozone Check	:No		
Ref. Male low	:0 U/I		
Ref. Male high	:11 U/I		
Ref. Female low	:0 U/I		
Ref. Female high	:11 U/I		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:30.0 ul
Rerun volume	:15.0 ul
R2 bottle	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Delay, min. time	:236, 283 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. Depletion	:0.200
Reagent Blank	:No
Factor	:853.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:30.0 ul
Rerun volume	:15.0 ul
Delay, min. time	:226, 196 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:0.200
Reagent Blank	:No
Factor	:853.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03



Analyzer: Selectra E
Test: ALT
Catalog # A7526

TEST PARAMETERS

Name	:ALT		
Abbr. Name	:ALT		
Mode	:Kinetic		
Wavelength	:340 nm		
Units	:U/I		
Decimals	:0	Ref. Ped. Low	:4 U/I
Low Conc.	:0 U/I	Ref. Ped. High	:36 U/I
High Conc.	:500 U/I	Control 1	:User Defined
Calibrator Name	:none	Control 2	:User Defined
Prozone Check	:No	Control 3	:User Defined
Ref. Male low	:4 U/I	Correlat. Factor	:1.000
Ref. Male high	:36 U/I	Correlat offset	:0.000
Ref. Female low	:4 U/I		
Ref. Female high	:36 U/I		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle	:5 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Delay, min. time	:50, 106 sec
Linearity Limit	:15.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.100
R. Abs. H. Limit	:2.500
Substr. Depletion	:0.500
Reagent Blank	:Yes
Factor	:-3376.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
Delay, min. time	:77, 159 sec
Linearity Limit	:15.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.800
R. Abs. H. Limit	:2.500
R. Abs. Deviation	:0.500
Reagent Blank	:Yes
Factor	:-3376.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E
 Test: ALP
 Catalog # A7516

TEST PARAMETERS

Name	: Alkaline Phosphatase	PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03
Abbr. Name	: ALP	
Mode	: Kinetic	
Wavelength	: 405 nm	Ref. Ped. Low : 35 U/l
Units	: U/l	Ref. Ped. High : 123 U/l
Decimals	: 0	Control 1 : User Defined
Low Conc.	: 0 U/l	Control 2 : User Defined
High Conc.	: 2000 U/l	Control 3 : User Defined
Calibrator Name	: none	Correlat. Factor : 1.000
Prozone Check	: No	Correlat offset : 0.000
Ref. Male low	: 35 U/l	
Ref. Male high	: 123 U/l	
Ref. Female low	: 35 U/l	
Ref. Female high	: 123 U/l	

DUAL MODE:

Sample Blank	: No
R1 bottle	: 25 ml
Normal volume	: 250 ul
Rerun volume	: 250 ul
Sample	
Normal volume	: 6.0 ul
Rerun volume	: 3.0 ul
R2 bottle	: 8 ml
Normal volume	: 50 ul
Rerun volume	: 50 ul
Predilution	: No
Slope Blank	: No
Delay, min. time	: 24, 106 sec
Linearity Limit	: 10.00%
Low absorbance	: 0.000
High absorbance	: 3.000
R. Abs. L. Limit	: 0.000
R. Abs. H. Limit	: 1.800
Substr. Depletion	: 3.000
Reagent Blank	: Yes
Factor	: 2720.00

MONO MODE:

Sample Blank	: No
R1 bottle	: 25 ml
Normal volume	: 250 ul
Rerun volume	: 250 ul
Sample	
Normal volume	: 5.0 ul
Rerun volume	: 2.5 ul
Delay, min. time	: 51, 117 sec
Linearity Limit	: 10.00%
Low absorbance	: 0.000
High absorbance	: 3.000
R. Abs. L. Limit	: 0.000
R. Abs. H. Limit	: 1.800
R. Abs. Deviation	: 3.000
Reagent Blank	: Yes
Factor	: 2720.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E
Test: Amylase
Catalog # A7553

TEST PARAMETERS

Name	:Amylase	PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03
Abbr. Name	:AMY	
Mode	:Kinetic	
Wavelength	:405 nm	
Units	:U/I	Ref. Ped. Low :25 U/I
Decimals	:0	Ref. Ped. High :125 U/I
Low Conc.	:0 U/I	Control 1 :User Defined
High Conc.	:2500 U/I	Control 2 :User Defined
Calibrator Name	:none	Control 3 :User Defined
Prozone Check	:No	Correlat. Factor :1.000
Ref. Male low	:25 U/I	Correlat offset :0.000
Ref. Male high	:125 U/I	
Ref. Female low	:25 U/I	
Ref. Female high	:125 U/I	

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:2.5 ul
R2 bottle	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Delay, min. time	:24, 132 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.600
Substr. Depletion	:3.000
Reagent Blank	:Yes
Factor	:3953.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:2.5 ul
Delay, min. time	:32, 136 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.600
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Factor	:3953.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E
Test: AST
Catalog # A7561

TEST PARAMETERS

Name	:AST	PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03
Abbr. Name	:AST	
Mode	:Kinetic	
Wavelength	:340 nm	
Units	:U/I	Ref. Ped. Low :5 U/I
Decimals	:0	Ref. Ped. High :34 U/I
Low Conc.	:0 U/I	Control 1 :User Defined
High Conc.	:500 U/I	Control 2 :User Defined
Calibrator Name	:none	Control 3 :User Defined
Prozone Check	:No	Correlat. Factor :1.000
Ref. Male low	:5 U/I	Correlat offset :0.000
Ref. Male high	:34 U/I	
Ref. Female low	:5 U/I	
Ref. Female high	:34 U/I	

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle	:5 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Delay, min. time	:50, 106 sec
Linearity Limit	:15.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.100
R. Abs. H. Limit	:2.500
Substr. Depletion	:0.500
Reagent Blank	:Yes
Factor	:-3376.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
Delay, min. time	:77, 159 sec
Linearity Limit	:15.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.800
R. Abs. H. Limit	:2.500
R. Abs. Deviation	:0.500
Reagent Blank	:Yes
Factor	:-3376.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E

Test: CK

Catalog # A7522

TEST PARAMETERS

Name	: Creatine Kinase		
Abbr. Name	: CK		
Mode	: Kinetic	Ref. Ped. Low	: 0 U/I
Wavelength	: 340 nm	Ref. Ped. High	: 160 U/I
Units	: U/I	Control 1	: User Defined
Decimals	: 0	Control 2	: User Defined
Low Conc.	: 0 U/I	Control 3	: User Defined
High Conc.	: 2000U/I	Correlat. Factor	: 1.000
Calibrator Name	: none	Correlat offset	: 0.000
Prozone Check	: No		
Ref. Male low	: 0 U/I		
Ref. Male high	: 160 U/I		
Ref. Female low	: 0 U/I		
Ref. Female high	: 130 U/I		

DUAL MODE:

Sample Blank	: No
R1 bottle	: 25 ml
Normal volume	: 250 ul
Rerun volume	: 250 ul
Sample	
Normal volume	: 15.0 ul
Rerun volume	: 7.5 ul
R2 bottle	: 8 ml
Normal volume	: 50 ul
Rerun volume	: 50 ul
Predilution	: No
Slope Blank	: No
Delay, min. time	: 103, 133 sec
Linearity Limit	: 10.00%
Low absorbance	: 0.000
High absorbance	: 3.000
R. Abs. L. Limit	: 0.000
R. Abs. H. Limit	: 0.700
Substr. Depletion	: 3.000
Reagent Blank	: Yes
Factor	: 3376.00

MONO MODE:

Sample Blank	: No
R1 bottle	: 25 ml
Normal volume	: 300 ul
Rerun volume	: 300 ul
Sample	
Normal volume	: 15.0 ul
Rerun volume	: 7.5 ul
Delay, min. time	: 129, 175 sec
Linearity Limit	: 10.00%
Low absorbance	: 0.000
High absorbance	: 3.000
R. Abs. L. Limit	: 0.000
R. Abs. H. Limit	: 0.700
R. Abs. Deviation	: 3.000
Reagent Blank	: Yes
Factor	: 3376.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03



Analyzer: Selectra E

Test: CK-MB

Catalog # C7562

TEST PARAMETERS

Name	:CK-MB		
Abbr. Name	:CKMB		
Mode	:Kinetic	Ref. Ped. Low	:0 U/I
Wavelength	:340 nm	Ref. Ped. High	:160 U/I
Units	:U/I	Control 1	:User Defined
Decimals	:0	Control 2	:User Defined
Low Conc.	:0 U/I	Control 3	:User Defined
High Conc.	:2000U/I	Correlat. Factor	:1.000
Calibrator Name	:none	Correlat offset	:0.000
Prozone Check	:No		
Ref. Male low	:0 U/I		
Ref. Male high	:160 U/I		
Ref. Female low	:0 U/I		
Ref. Female high	:160 U/I		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:30.0 ul
Rerun volume	:30.0 ul
R2 bottle	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Delay, min. time	:307, 212 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.700
Substr. Depletion	:3.000
Reagent Blank	:Yes
Factor	:3537.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:30.0 ul
Rerun volume	:30.0 ul
Delay, min. time	:285, 137 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:0.700
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Factor	:3537.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03



Analyzer: Selectra E
Test: GGT
Catalog # G7571

TEST PARAMETERS

Name	:GGT		
Abbr. Name	:GGT		
Mode	:Kinetic	Ref. Ped. Low	:8 U/l
Wavelength	:405 nm	Ref. Ped. High	:54 U/l
Units	:U/l	Control 1	:User Defined
Decimals	:0	Control 2	:User Defined
Low Conc.	:0 U/l	Control 3	:User Defined
High Conc.	:1000 U/l	Correlat. Factor	:1.000
Calibrator Name	:none	Correlat offset	:0.000
Prozone Check	:No		
Ref. Male low	:8 U/l		
Ref. Male high	:54 U/l		
Ref. Female low	:8 U/l		
Ref. Female high	:54 U/l		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:7.5 ul
R2 bottle	:8 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Delay, min. time	:24, 132 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. Depletion	:3.000
Reagent Blank	:Yes
Factor	:2188.00

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:7.5 ul
Delay, min. time	:32, 194 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Factor	:2188.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E
 Test: Lactate
 Catalog # L7596

TEST PARAMETERS

Name	:Lactate	Reagent Blank	:Yes
Dehydrogenase		Factor	:6190.0
Abbr. Name	:LD		
Mode	:Kinetic		
Wavelength	:340 nm	Ref. Ped. Low	:80 U/l
Units	:U/l	Ref. Ped. High	:285 U/l
Decimals	:0	Control 1	:User Defined
Low Conc.	:0 U/l	Control 2	:User Defined
High Conc.	:1500 U/l	Control 3	:User Defined
Calibrator Name	:none	Correlat. Factor	:1.000
Prozone Check	:No	Correlat offset	:0.000
Ref. Male low	:80 U/l		
Ref. Male high	:285 U/l		
Ref. Female low	:80 U/l		
Ref. Female high	:285 U/l		

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:8.0 ul
Rerun volume	:4.0 ul
R2 bottle	:8 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Delay, min. time	:24, 79 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.600
Substr. Depletion	:3.000

MONO MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:8.0 ul
Rerun volume	:4.0 ul
Delay, min. time	:32, 78 sec
Linearity Limit	:10.00%
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.600
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Factor	:6190.00

Notes: It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50. Rev: 1-03

Pointe Scientific Test Parameters - Selectra E

Creatinine (CREA): C7539

TEST PARAMETERS

Name	:Creatinine	Prozone Check	:No
Abbr. Name	:CREA	Ref. Male low	:0.4 mg/dl
Mode	:Twopoint	Ref. Male high	:1.4 mg/dl
Wavelength	:505 nm	Ref. Female low	:0.4 mg/dl
Units	:mg/dl	Ref. Female high	:1.4 mg/dl
Decimals	:1	Ref. Ped. Low	:0.4 mg/dl
Low Conc.	:0.0 mg/dl	Ref. Ped. High	:1.4 mg/dl
High Conc.	:25.0 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Point one, two	:50, 130 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.400
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
Point one, two	:51, 129 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:0.400
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Fructosamine (FRUC): F7546

TEST PARAMETERS

Name	:Fructosamine	Prozone Check	:No
Abbr. Name	:FRUC	Ref. Male low	:1.61 mmol/l
Mode	:Twopoint	Ref. Male high	:2.68 mmol/l
Wavelength	:546 nm	Ref. Female low	:1.61 mmol/l
Units	:mmol/l	Ref. Female high	:2.68 mmol/l
Decimals	:2	Ref. Ped. Low	:1.61 mmol/l
Low Conc.	:0.00 mmol/l	Ref. Ped. High	:2.68 mmol/l
High Conc.	:8.50 mmol/l	Control 1	:User Defined
Calibrator Name	:FrucCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Point one, two	:307, 519 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
Point one, two	:285, 422 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Lipase (LIP): L7503

TEST PARAMETERS

Name	:Lipase	Prozone Check	:No
Abbr. Name	:LIP	Ref. Male low	:0 U/I
Mode	:Twopoint	Ref. Male high	:62 U/I
Wavelength	:546 nm	Ref. Female low	:0 U/I
Units	:U/I	Ref. Female high	:62 U/I
Decimals	:0	Ref. Ped. Low	:0 U/I
Low Conc.	:0 U/I	Ref. Ped. High	:62 U/I
High Conc.	:600 U/I	Control 1	:User Defined
Calibrator Name	:LIP STD	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:295 ul
Rerun volume	:295 ul
Sample	
Normal volume	:5.0 ul
Rerun volume	:5.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:100 ul
Rerun volume	:100 ul
Predilution	:No
Slope Blank	:No
Point one, two	:130, 236 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:-0.100
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:2.000
Reagent Blank	:No
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

MONO MODE:

Sample Blank	:
R1 bottle (ml)	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

Urea Nitrogen (BUN): B7552

TEST PARAMETERS

Name	:Urea Nitrogen	Prozone Check	:No
Abbr. Name	:BUN	Ref. Male low	:7 mg/dl
Mode	:Twopoint	Ref. Male high	:18 mg/dl
Wavelength	:340 nm	Ref. Female low	:7 mg/dl
Units	:mg/dl	Ref. Female high	:18 mg/dl
Decimals	:0	Ref. Ped. Low	:7 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. High	:18 mg/dl
High Conc.	:150 mg/dl	Control 1	:User Defined
Calibrator Name	:ChemCal	Control 2	:User Defined
Repeat	:1	Control 3	:User Defined
Number	:1	Correlat. Factor	:1.000
Concentration	:User Defined	Correlat offset	:0.000
Interval	:7 days		
Cut-off	:No		

DUAL MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
R2 bottle (ml)	:8 ml
Normal volume	:
Rerun volume	:
Predilution	:No
Slope Blank	:No
Point one, two	:24, 77 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:1.000
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

MONO MODE:

Sample Blank	:No
R1 bottle (ml)	:25 ml
Normal volume	:300 ul
Rerun volume	:300 ul
Sample	
Normal volume	:3.0 ul
Rerun volume	:3.0 ul
Point one, two	:32, 90 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:1.000
R. Abs. H. Limit	:3.000
R. Abs. Deviation	:3.000
Reagent Blank	:Yes
Cal. Low Limit	:0.000 dAbs
Cal. High Limit	:0.000 dAbs
Factor	:1.00

Notes:

Pointe Scientific Test Parameters - Selectra E

Apolipoprotein A1 (APOA): A7544

TEST PARAMETERS

Name :Apolipoprotein A1
Abbr. Name :APOA
Mode :Twopoint
Wavelength :340 nm
Units :mg/dl
Decimals :0
Low Conc. :0 mg/dl
High Conc. :230 mg/dl
Calibrator Name :APO A1 CAL
Repeat :1
Number :6
Concentration :User Defined
Interval :7 days
Cut-off :1.00%

Notes: Mono Mode not applicable

Rev. 3/00
Prozone Check :No
Ref. Male low :107 mg/dl
Ref. Male high :187 mg/dl
Ref. Female low :115 mg/dl
Ref. Female high :206 mg/dl
Ref. Ped. Low :107 mg/dl
Ref. Ped. High :206 mg/dl
Control 1 :User Defined
Control 2 :User Defined
Control 3 :User Defined
Correlat. Factor :1.000
Correlat offset :0.000

DUAL MODE:

Sample Blank :No
R1 bottle (ml) :25 ml
Normal volume :250 ul
Rerun volume :250 ul
Sample
Normal volume :4.0 ul
Rerun volume :4.0 ul
R2 bottle (ml) :8 ml
Normal volume :135 ul
Rerun volume :135 ul
Predilution :No
Slope Blank :No
Point one, two :-3, 236 sec
Low absorbance :0.000
High absorbance :3.000
R. Abs. L. Limit :0.000
R. Abs. H. Limit :3.000
R. Abs. Deviation :2.300
Reagent Blank :No
Cal. Low Limit :0.000 dAbs
Cal. High Limit :0.000 dAbs
Factor :N/A

MONO MODE:

Sample Blank :
R1 bottle (ml) :
Normal volume :
Rerun volume :
Sample
Normal volume :
Rerun volume :
Point one, two :
Low absorbance :
High absorbance :
R. Abs. L. Limit :
R. Abs. H. Limit :
R. Abs. Deviation :
Reagent Blank :
Cal. Low Limit :0.000 dAbs
Cal. High Limit :0.000 dAbs
Factor :

Pointe Scientific Test Parameters - Selectra E

Apolipoprotein B (APOB): A7588

TEST PARAMETERS

Name	:Apolipoprotein B	Max. inaccuracy	:1.00%
Abbr. Name	:APOB	Prozone Check	:No
Mode	:Twopoint	Ref. Male low	:63 mg/dl
Wavelength	:340 nm	Ref. Male high	:152 mg/dl
Units	:mg/dl	Ref. Female low	:60 mg/dl
Decimals	:0	Ref. Female high	:150 mg/dl
Low Conc.	:0 mg/dl	Ref. Ped. Low	:60 mg/dl
High Conc.	:215 mg/dl	Ref. Ped. High	:152 mg/dl
Calibrator Name	:APO B CAL	Control 1	:User Defined
Repeat	:1	Control 2	:User Defined
Number	:6	Control 3	:User Defined
Interval	:7 days	Correlat. Factor	:1.000
Cut-off	:No	Correlat offset	:0.000

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:4.0 ul
Rerun volume	:4.0 ul
R2 bottle	:8 ml
Normal volume	:135 ul
Rerun volume	:135 ul
Predilution	:No
Slope Blank	:No
Point one, two	:-3, 236 sec
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:1.200
Reagent Blank	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

ASO (ASO): A7566

TEST PARAMETERS

Name	:ASO	Max. inaccuracy	:1.00%
Abbr. Name	:ASO	Prozone Check	:No
Mode	:Twopoint	Ref. Male low	:0 U/ml
Wavelength	:340 nm	Ref. Male high	:156 U/ml
Units	:U/ml	Ref. Female low	:0 U/ml
Decimals	:0	Ref. Female high	:156 U/ml
Low Conc.	:0 U/ml	Ref. Ped. Low	:0 U/ml
High Conc.	:1000 U/ml	Ref. Ped. High	:156 U/ml
Calibrator Name	:ASO CAL	Control 1	:User Defined
Repeat	:1	Control 2	:User Defined
Number	:6	Control 3	:User Defined
Interval	:7 days	Correlat. Factor	:1.000
Cut-off	:No	Correlat offset	:0.000

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle	:8 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Point one, two	:-3, 236 sec
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:0.500
Reagent Blank	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

C-Reactive Protein HS (CRP-HS): C7564

TEST PARAMETERS

Name	:C-Reactive Protein HS	Rev. 3/00	
Abbr. Name	:CRP-HS	Max. inaccuracy	:1.00%
Mode	:Twopoint	Prozone Check	:No
Wavelength	:546nm	Ref. Male low	:0.0 mg/dl
Units	:mg/dl	Ref. Male high	:0.6 mg/dl
Decimals	:3	Ref. Female low	:0.0 mg/dl
Low Conc.	:0.005mg/dl	Ref. Female high	:0.6 mg/dl
High Conc.	:1.000mg/dl	Ref. Ped. Low	:0.0 mg/dl
Calibrator Name	:CRP CAL	Ref. Ped. High	:0.6 mg/dl
Repeat	:1	Control 1	:User Defined
Number	:6	Control 2	:User Defined
Interval	:7 days	Control 3	:User Defined
Cut-off	:No	Correlat. Factor	:1.000
		Correlat offset	:0.000

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:150 ul
Rerun volume	:150 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle	:5 ml
Normal volume	:100 ul
Rerun volume	:100 ul
Predilution	:No
Slope Blank	:No
Point one, two	:6, 236 sec
Low absorbance	:-0.100
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:2.000
Reagent Blank	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

Microalbumin (UALB): M7562

TEST PARAMETERS

Name	:Microalbumin	Max. inaccuracy	:1.00%
Abbr. Name	:UALB	Prozone Check	:No
Mode	:Twopoint	Ref. Male low	:0.0 mg/dl (24h)
Wavelength	:340 nm	Ref. Male high	:2.1 mg/dl (24h)
Units	:mg/dl	Ref. Female low	:0.0 mg/dl (24h)
Decimals	:1	Ref. Female high	:2.1 mg/dl (24h)
Low Conc.	:0.0 mg/dl	Ref. Ped. Low	:0.0 mg/dl (24h)
High Conc.	:30.0 mg/dl	Ref. Ped. High	:2.1 mg/dl (24h)
Calibrator Name	:UALB CAL	Control 1	:User Defined
Repeat	:1	Control 2	:User Defined
Number	:6	Control 3	:User Defined
Interval	:7 days	Correlat. Factor	:1.000
Cut-off	:No	Correlat offset	:0.000

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:285 ul
Rerun volume	:285 ul
Sample	
Normal volume	:12.0 ul
Rerun volume	:12.0 ul
R2 bottle	:8 ml
Normal volume	:95 ul
Rerun volume	:95 ul
Predilution	:No
Slope Blank	:No
Point one, two	:6, 236 sec
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:1.400
Reagent Blank	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

Rheumatoid Factor (RF): R7568

TEST PARAMETERS

Name	:Rheumatoid Factor	Max. inaccuracy	:1.00%
Abbr. Name	:RF	Prozone Check	:No
Mode	:Twopoint	Ref. Male low	:0 U/ml
Wavelength	:340 nm	Ref. Male high	:10 U/ml
Units	:U/ml	Ref. Female low	:0 U/ml
Decimals	:0	Ref. Female high	:10 U/ml
Low Conc.	:0 U/ml	Ref. Ped. Low	:0 U/ml
High Conc.	:300 U/ml	Ref. Ped. High	:10 U/ml
Calibrator Name	:RF CAL	Control 1	:User Defined
Repeat	:1	Control 2	:User Defined
Number	:4	Control 3	:User Defined
Interval	:7 days	Correlat. Factor	:1.000
Cut-off	:No	Correlat offset	:0.000

DUAL MODE:

Sample Blank	:No
R1 bottle	:25 ml
Normal volume	:250 ul
Rerun volume	:250 ul
Sample	
Normal volume	:15.0 ul
Rerun volume	:15.0 ul
R2 bottle	:8 ml
Normal volume	:50 ul
Rerun volume	:50 ul
Predilution	:No
Slope Blank	:No
Point one, two	:-3, 236 sec
Low absorbance	:0.000
High absorbance	:3.000
R. Abs. L. Limit	:0.000
R. Abs. H. Limit	:3.000
Substr. depletion	:0.500
Reagent Blank	:No

MONO MODE:

Sample Blank	:
R1 bottle	:
Normal volume	:
Rerun volume	:
Sample	
Normal volume	:
Rerun volume	:
Point one, two	:
Low absorbance	:
High absorbance	:
R. Abs. L. Limit	:
R. Abs. H. Limit	:
R. Abs. Deviation	:
Reagent Blank	:

Notes: Mono Mode not applicable

Pointe Scientific Test Parameters - Selectra E

UIBC (UIBC): 17506

TEST PARAMETERS

Name :UIBC
Abbr. Name :UIBC
Mode :
Wavelength :546 nm
Units :ug/dl
Decimals :0
Low Conc. :
High Conc. :
Calibrator Name :ChemCal
Repeat :1
Number :1
Interval :7 days
Cut-off :No

DUAL MODE:

Sample Blank :No
R1 bottle :25 ml
Normal volume :
Rerun volume :
Sample
Normal volume :
Rerun volume :
R2 bottle :8 ml
Normal volume :
Rerun volume :
Predilution :No
Slope Blank :
Point one, two :
Low absorbance :0.000
High absorbance :3.000
R. Abs. L. Limit :
R. Abs. H. Limit :1.000
Substr. depletion :
Reagent Blank :Yes

Notes:

Rev. 3/00
Max. inaccuracy :No
Prozone Check :No
Ref. Male low :
Ref. Male high :
Ref. Female low :
Ref. Female high :
Ref. Ped. Low :
Ref. Ped. High :
Control 1 :User Defined
Control 2 :User Defined
Control 3 :User Defined
Correlat. Factor :1.000
Correlat offset :0.000

MONO MODE:

Sample Blank :No
R1 bottle :25 ml
Normal volume :
Rerun volume :
Sample
Normal volume :
Rerun volume :
Point one, two :
Low absorbance :
High absorbance :
R. Abs. L. Limit :
R. Abs. H. Limit :
R. Abs. Deviation :
Reagent Blank :Yes

Pointe Scientific, Inc.

Instrument Application

Analyzer: Selectra E
Test: Lactate Dehydrogenase
Catalog # L7572

TEST PARAMETERS

Name : Lactate Dehydrogenase
 Abbr. Name : LD
 Mode : Kinetic
 Wavelength : 340 nm
 Units : U/l
 Decimals : 0
 Low Conc. : 0 U/l
 High Conc. : 1500 U/l
 Calibrator Name : none
 Repeat :
 Number :
 Concentration : User Defined
 Interval : N/A
 Cut-off : N/A

Prozone Check : No
 Ref. Male low : 80 U/l
 Ref. Male high : 285 U/l
 Ref. Female low : 80 U/l
 Ref. Female high : 285 U/l
 Ref. Ped. Low : 80 U/l
 Ref. Ped. High : 285 U/l
 Control 1 : User Defined
 Control 2 : User Defined
 Control 3 : User Defined
 Correlat. Factor : 1.000
 Correlat offset : 0.000

DUAL MODE:

Sample Blank : No
 R1 bottle (ml) : 25 ml
 Normal volume : 250 ul
 Rerun volume : 250 ul
 Sample
 Normal volume : 8.0 ul
 Rerun volume : 4.0 ul
 R2 bottle (ml) : 8 ml
 Normal volume : 50 ul
 Rerun volume : 50 ul
 Predilution : No
 Incubation time : 24, 79 sec
 Low absorbance : 0.000
 High absorbance : 3.000
 R. Abs. L. Limit : -0.100
 R. Abs. H. Limit : 0.600
 Reagent Blank : Yes
 Cal. Low Limit : 0.000 Abs
 Cal. High Limit : 0.000 Abs
 Factor : 6190.00

MONO MODE:

Sample Blank : No
 R1 bottle (ml) : 25 ml
 Normal volume : 300 ul
 Rerun volume : 300 ul
 Sample
 Normal volume : 8.0 ul
 Rerun volume : 4.0 ul
 Incubation time : 32, 78 sec
 Low absorbance : 0.000
 High absorbance : 3.000
 R. Abs. L. Limit : -0.100
 R. Abs. H. Limit : 0.600
 Reagent Blank : Yes
 Cal. Low Limit : 0.000 Abs
 Cal. High Limit : 0.000 Abs
 Factor : 6190.00

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.
 Rev 1-03