

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** ASO  
**Catalog # :** A7566

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: LIN. INTER  
 Reagent Blank: NO  
 Cleaner: NO (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/ML

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 20.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 250 UL  
 Start Reag. 1 Cycle: 10  
 Volume: 50 UL  
 Diluent Name: H2O  
 Volume: 10 UL

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 500 U/L  
 Normal Range Low: \* U/ML  
 Normal Range High: \* U/ML

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 9  
 Reading Last: 25  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: NO  
 Reag. Range High: NO  
 Blank Range Low: NO  
 Blank Range High: NO  
 Factor: -  
 Standard Pos: \*  
 STD-1: 0.0 (a)  
 STD-2-6:

Replicate: Single  
 Deviation: NO

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Use Saline as 0.0 std.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 11-03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test:  $\beta$ - Hydroxybutyrate  
 Catalog # : H7587

**Reagent preparation:** Reagents provided as ready to use liquids.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S -SR1  
 Calibration Mode: Slope Avg.  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: NO  
 Wavelength: 500nm  
 Decimal Position: 2  
 Unit: mmol/L

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 uL  
 Diluent Name: H2O  
 Volume: 20.0 uL

Reagent Cycle: 1  
 Volume: 105 uL  
 Start Reag. 1 Cycle: 2  
 Volume: 18 uL  
 Diluent: H2O  
 Volume: 60.0 uL

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0  
 Test Range High: 4.5  
 Normal Range Low: 0.02  
 Normal Range High: 0.27

Number of Steps: 1

Calc. Step A: Endpoint  
 Reading First: 1  
 Reading Last: 12  
 Reaction Limit: NO  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: NO  
 Reag. Range High: NO  
 Blank Range Low: NO  
 Blank Range High: NO  
 Factor: -

Standard Pos: \*  
 STD-1: \*  
 STD-2:  
 STD-3:

Replicate: Single  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED  
 Rev: 3-03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Acid Phosphatase  
 Catalog # : A7503

Add 8.5 ml deionized water to the 10.0 ml vials.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil  
 Cleaner: Before (2)

Wavelength: 405nm  
 Decimal Position: 1  
 Unit: U/L

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 20.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 150 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On  
 Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: 0.0 U/L  
 Test Range High: 35.0 U/L  
 Normal Range Low: 0.0 U/L  
 Normal Range High: 9.0 U/L

Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 13  
 Reading Last: 25  
 Reaction Limit: 0.400  
 Point: T1

### CALIBRATION

Calib Interval: On Request  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: 1421

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 1/03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Albumin**  
**Catalog # : A7502**

Reagent provided ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 600nm (5)  
 Decimal Position: 1  
 Unit: G/DL (11)

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.050  
 Reag. Range High: 0.450  
 Blank Range Low: -0.010  
 Blank Range High: 0.375  
 Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 2.0 UL  
 Diluent Name: H2O  
 Volume: 28.0 UL

Reagent Cycle: 1  
 Volume: 370 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

\* USER DEFINED

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

Convers. Factor: 1.00000  
 Offset: 0.00000

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Test Range Low: 0.0 G/DL  
 Test Range High: 6.0 G/DL  
 Normal Range Low: 3.5 G/DL  
 Normal Range High: 5.3 G/DL

Rev. 1-03

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: CB  
 Reading Last: T2  
 Reaction Limit: -  
 Point: -

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Alcohol**  
**Catalog # : A7504**

Prepare reagent according to package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Std Lin or Cal  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: After  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mg/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 10.0 UL  
 Diluent Name: H2O  
 Volume: 40.0 UL

Reagent Cycle: 1  
 Volume: 260 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit:  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0  
 Test Range High: 400  
 Normal Range Low: 0

Normal Range High: User Defined

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 17  
 Reaction Limit: No

### CALIBRATION

Calib Interval: Each Run  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: 1.500  
 Blank Range Low: No  
 Blank Range High: 1.500  
 Factor: -

Standard Pos: -  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos: \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos: \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Alcohol Standard A7504-STD recommended for calibration.

\* USER DEFINED  
 Rev: 1-03

It is recommended that two levels of control material be assayed daily. Reorder PSI Controls Cat.# A7504-CTL.

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Alkaline Phosphatase (Liquid)  
**Catalog # :** A7516

Single working reagent is prepared by mixing 4 parts R1 with 1 part R2.

	<u>Single Reagent</u>	<u>2Part App.</u>
<b><u>GENERAL</u></b>		
Measurement Mode:	Absorb	
Reaction Mode:	R-S-(1)	R-S-SR1(3)
Calibration Mode:	Factor (1)	
Reagent Blank:	Reag/Dil (2)	
Cleaner:	No (1)	
Wavelength:	405nm (2)	
Decimal Position:	0	
Unit:	U/L (21)	

<b><u>ANALYSIS</u></b>		
Post Dil. Factor:	No	
Post Conc. Factor:	No	
Sample Cycle:	2	2
Volume:	4.0 UL	4.0
Diluent Name:	H2O	H2O
Volume:	30.0 UL	20
Reagent Cycle:	1	1
Volume:	170 UL	135
Start Reag. 1 Cycle:	-	1
Volume:	-	35
Diluent Name:	-	H2O
Volume:	-	10.0

<b><u>CALCULATION</u></b>		
Sample Limit:	No	
Point:	-	
Reac. Direction:	Increase (1)	
Check:	Off	

Convers. Factor:	1.00000
Offset:	0.00000

Test Range Low:	0 U/L
Test Range High:	1500 U/L
Normal Range Low:	35 U/L
Normal Range High:	123 U/L

Number of Steps:	1
------------------	---

Calc. Step A:	Kinsearch (3)
Reading First:	5
Reading Last:	10

Reaction Limit:	No
Point:	-

<b><u>CALIBRATION</u></b>	
Calib Interval:	On Request (3)
Blank	
Reag. Range Low:	0.1200
Reag. Range High:	0.5000
Blank Range Low:	-0.0050
Blank Range High:	0.0050

Factor:	4521
Standard Pos:	-
STD-1:	-
STD-2:	
STD-3:	

Replicate:	Dupl (2)
Deviation:	10%

<b><u>Control</u></b>	
CS1 Pos	* Low: (User Defined) Assign: High:
CS2 Pos	* Low: (User Defined) Assign: High:

CS3 Pos:	No
----------	----

\*USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 1-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Alkaline Phosphatase**  
**Catalog # : A7505**

Add 12ml and 40ml to 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 405nm (2)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.3500  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 1500 U/L  
 Normal Range Low: 35 U/L  
 Normal Range High: 123 U/L

Number of Steps: 1  
 Calc. Step A: Kinsearch (3)  
 Reading First: 2  
 Reading Last: 6  
 Reaction Limit: 0.530  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0500  
 Reag. Range High: 0.8000  
 Blank Range Low: -.0040  
 Blank Range High: 0.0080

Factor: 2844

Standard Pos: -

STD-1: -

STD-2:

STD-3:

Replicate: Dupl (2)

Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)

Assign:

High:

CS2 Pos \* Low: (User Defined)

Assign:

High:

CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 1-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: ALT (SGPT) Liquid**  
**Catalog # : A7526**

Single working reagent is prepared by mixing 4 parts R1 with 1 part R2.

Single reagent      2Part app.

**GENERAL**

Measurement Mode:	Absorb	
Reaction Mode:	R-S-(1)	R-S-SR1 (3)
Calibration Mode:	Factor (1)	
Reagent Blank:	Reag/Dil (2)	
Cleaner:	Before (2)	
Wavelength:	340nm (1)	
Decimal Position:	0	
Unit:	U/L (21)	

**ANALYSIS**

Post Dil. Factor:	No	
Post Conc. Factor:	No	
Sample Cycle:	2	2
Volume:	12.0 UL	12.0
Diluent Name:	H2O	H2O
Volume:	25.0 UL	10.0
Reagent Cycle:	1	1
Volume:	155 UL	125.0
Start Reag. 1 Cycle:	-	1
Volume:	-	30.0
Diluent Name:	-	H2O
Volume:	-	15.0

**CALCULATION**

Sample Limit:	No
Point:	-
Reac. Direction:	Decrease (2)
Check:	On
Convers. Factor:	1.00000
Offset:	0.00000
Test Range Low:	0 U/L
Test Range High:	600 U/L
Normal Range Low:	4 U/L
Normal Range High:	36 U/L
Number of Steps:	1
Calc. Step A:	Kinsearch (3)
Reading First:	6
Reading Last:	12
Reaction Limit:	No
Point:	-

**CALIBRATION**

Calib Interval:	On Request (3)
Blank	
Reag. Range Low:	0.4200
Reag. Range High:	2.0000
Blank Range Low:	-0.0050
Blank Range High:	0.0050
Factor:	4788
Standard Pos:	-
STD-1:	-
STD-2:	
STD-3:	

Replicate:	Dupl (2)
Deviation:	10%

**Control**

CS1 Pos	* Low: (User Defined)
	Assign:
	High:
CS2 Pos	* Low: (User Defined)
	Assign:
	High:
CS3 Pos:	No

\*USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 5/03



# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** ALT (SGPT)  
**Catalog # :** A7525

Add 12ml and 40ml DiH2O to 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 13.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.7000  
 Point: T1  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 500 U/L  
 Normal Range Low: 0 U/L  
 Normal Range High: 38 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 11  
 Reaction Limit: .1450  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.5400  
 Reag. Range High: 1.6000  
 Blank Range Low: -0.0040  
 Blank Range High: 0.0040  
 Factor: 3463

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 1-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Ammonia  
**Catalog # :** A7553

Add 5mls ammonia free DH2O to the 6.5ml substrate rgt. (R1)  
 Add 2mls of ammonia free DH2O to the enzyme reagent vial. (SR1)

**GENERAL**

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Std. Lin.  
 Reagent Blank: Reag/Dil  
 Cleaner: No  
 Wavelength: 340nm  
 Decimal Position: 0  
 Unit: umol/L (ug/dl)

**CALIBRATION**

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No  
 Factor: -

**ANALYSIS**

Post Dil. Factor: No  
 Post Conc. Factor: No

Standard Pos:  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Sample Cycle: 1  
 Volume: 40.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Replicate: Dupl (2)  
 Deviation: 10%

Reagent Cycle: 1  
 Volume: 150 UL  
 Start Reag. 1 Cycle: 10  
 Volume: 8.0 UL  
 Diluent Name: H2O  
 Volume: 4.0 UL

**Control**  
 CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

**CALCULATION**

Sample Limit: No  
 Point: -  
 Reac. Direction: Decrease  
 Check: On

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Controls Cat.# A7504-CTL

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 600 U/L  
 Normal Range Low: \*  
 Normal Range High: \*  
 Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 9  
 Reading Last: 21  
 Reaction Limit: -  
 Point: -

Rev: 2/03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Amylase  
**Catalog # :** A7564

Reagent provided as a ready to use liquid.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 405 nm (2)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.2500  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 2000 U/L  
 Normal Range Low: 25 U/L  
 Normal Range High: 125 U/L

Number of Steps: 1  
 Calc. Step A: Kinsearch (3)  
 Reading First: 2  
 Reading Last: 6  
 Reaction Limit: NO  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0070  
 Reag. Range High: 0.5400  
 Blank Range Low: -0.0060  
 Blank Range High: 0.0050  
 Factor: 6804

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 2/02

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Apolipoprotein A1  
**Catalog # :** A7544

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Std Nonlin  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No  
 Calibrator Cup Pos: 1 -  
 Cal-1: \*  
 Cal 2-6: User Defined

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 2.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Replicate: Single  
 Deviation: NO

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: 2  
 Volume: 75 UL  
 Diluent: 10.0 UL

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: ON

Values exceeding the highest calibrator must be diluted with saline.

\* User defined. Saline should be used as the 0.0 standard

Convers. Factor: 1.00000  
 Offset: 0.00000

It is recommended that two levels of control material be assayed daily. Reorder PSI Lipid Controls Cat.# L7580-18

Test Range Low: 0.0 mg/dl  
 Test Range High: 200 mg/dl  
 Normal Range Low: User Defined  
 Normal Range High: User Defined

Rev: 11-03

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 12  
 Reaction Limit: -  
 Point: -

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Apolipoprotein B**  
**Catalog # : A7588**

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Std Nonlin  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No  
 Calibrator Cup Pos: 1 -  
 Cal-1: \*  
 Cal 2-6: User Defined

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Replicate: Single  
 Deviation: NO

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: 2  
 Volume: 75 UL  
 Diluent: 10.0 UL

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

### CALCULATION

Sample Limit: No  
 Point: -  
 Rec. Direction: Increase (1)  
 Check: On

Values exceeding the highest calibrator must be diluted with saline.

\* User defined. Saline should be used as the 0.0 standard.

Convers. Factor: 1.00000  
 Offset: 0.00000

It is recommended that two levels of control material be assayed daily. Reorder PSI Lipid Controls Cat.# L7580-18

Test Range Low: 0.0 mg/dl  
 Test Range High: 200 mg/dl  
 Normal Range Low: User Defined  
 Normal Range High: User Defined

Rev. 11-03

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 14  
 Reaction Limit: -  
 Point: -

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** ASO  
**Catalog # :** A7566

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: LIN. INTER  
 Reagent Blank: NO  
 Cleaner: NO (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/ML

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 20.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 250 UL  
 Start Reag. 1 Cycle: 10  
 Volume: 50 UL  
 Diluent Name: H2O  
 Volume: 10 UL

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 500 U/L  
 Normal Range Low: \* U/ML  
 Normal Range High: \* U/ML

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 9  
 Reading Last: 25  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: NO  
 Reag. Range High: NO  
 Blank Range Low: NO  
 Blank Range High: NO  
 Factor: -  
 Standard Pos: \*  
 STD-1: 0.0 (a)  
 STD-2-6:

Replicate: Single  
 Deviation: NO

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Use Saline as 0.0 std.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 11-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** AST (SGOT) Liquid  
**Catalog # :** A7561

Single working reagent is prepared by mixing 4 parts R1 with 1 part R2.

Single reagent    2Part app.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)    R-S-SR1 (3)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2                    2  
 Volume: 12.0 UL            12.0 UL  
 Diluent Name: H2O            H2O  
 Volume: 25.0 UL            10.0 UL

Reagent Cycle: 1                    1  
 Volume: 155 UL            125 UL  
 Start Reag 1 Cycle: -            1  
 Volume: -                    30.0 UL  
 Diluent Name: -                    H2O  
 Diluent: -                    15.0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 600 U/L  
 Normal Range Low: 5 U/L  
 Normal Range High: 34 U/L

Number of Steps: 1

Calc. Step A: Kinetic  
 Reading First: 6  
 Reading Last: 12  
 Reaction Limit: No  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.4200  
 Reag. Range High: 2.0000  
 Blank Range Low: -0.0050  
 Blank Range High: 0.0050

Factor: 4788

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos            \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos            \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos:            No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 5-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: AST (SGOT)**  
**Catalog # : A7560**

Add 12ml and 40ml to 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 13 UL  
 Diluent Name: H2O  
 Volume: 30 UL

Reagent Cycle: 1  
 Volume: 125 UL

### CALCULATION

Sample Limit: 0.7000  
 Point: T1  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 500 U/L  
 Normal Range Low: 0 U/L  
 Normal Range High: 40 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 11  
 Reaction Limit: .1450  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.5400  
 Reag. Range High: 1.6000  
 Blank Range Low: -0.0040  
 Blank Range High: 0.0040

Factor: 3463

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03



# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** auto HDL  
**Catalog # :** H7545

Reagents are supplied ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No  
  
 Wavelength: 600nm (5)  
 Decimal Position: 0  
 Unit: mg/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No  
  
 Sample Cycle: 1  
 Volume: 3.0 U/L  
 Diluent Name: H2O  
 Volume: 10.0 UL  
  
 Reagent Cycle: 1  
 Volume: 240 UL  
 Start Reag 1 Cycle: 12  
 Volume: 80.0 UL  
 Diluent Name: H2O  
 Diluent: 5.0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off  
  
 Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: No  
 Test Range High: No  
 Normal Range Low: 30 mg/dl  
 Normal Range High: 85 mg/dl  
  
 Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 11  
 Reading Last: 24  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No  
  
 Factor: -  
  
 Standard Pos: 1  
     STD-1: User Defined  
     STD-2:  
     STD-3:

Replicate: Dupl (2)  
 Deviation: No

### Control

CS1 Pos      \* Low: (User Defined)  
                   Assign:  
                   High:  
 CS2 Pos      \*Low: (User Defined)  
                   Assign:  
                   High:  
 CS3 Pos:      No

autoHDL/LDL calibrator recommended for calibration  
 Catalog# H7545-CAL

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 10/02

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** auto LDL  
**Catalog # :** L7574

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 550nm (4)  
 Decimal Position: 0  
 Unit: mg/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 2.4 UL  
 Diluent Name: H2O  
 Volume: 10 UL

Reagent Cycle: 1  
 Volume: 240 UL  
 Start Reag 1 Cycle: 12  
 Volume: 80.0 UL  
 Diluent Name: H2O  
 Diluent: 5.0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Rec. Direction: Increase (1)  
 Check: ON

Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: 0  
 Test Range High: 700  
 Normal Range Low: 20 mg/dl  
 Normal Range High: 130 mg/dl

Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 11  
 Reading Last: 24  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.100  
 Reag. Range High: 2.000  
 Blank Range Low: -0.1000  
 Blank Range High: 1.5000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Enter calibrator value. AutoHDL/LDL calibrator catalog# H7545-CAL is recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Bilirubin, Direct**  
**Catalog # : B7538**

Reagents are ready to use. NOTE: This is a two (2) reagent system used on (rack3). Direct Bilirubin reagent is placed in the larger container. Nitrite reagent in the smaller.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1 (3)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 30.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag 1 Cycle: 2  
 Volume: 4.0 UL  
 Diluent Name: H2O  
 Diluent: 20.0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 0.0 mg/dl  
 Normal Range High: 0.5 mg/dl

Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 2  
 Reading Last: 8

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -

Standard Pos: User Defined  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Bilirubin, Total**  
**Catalog # : B7576**

Prepare reagent according to package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor:

Sample Cycle: 2  
 Volume: 7.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 175 UL  
 Start Reag 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No

Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 0.2 mg/dl  
 Normal Range High: 1.2 mg/dl

Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 13

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0800  
 Reag. Range High: 0.1500  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor:

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

\* USER DEFINED

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** BUN, Liquid  
**Catalog # :** B7552

Reagents are prepared by mixing 5 parts R1 with 1 part R2.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
  
 Wavelength: 340nm (1)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 140.0 mg/dl  
 Normal Range Low: 7.0 mg/dl  
 Normal Range High: 18.0 mg/dl

Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 2  
 Reading Last: 6  
 Reaction Limit: No  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.5000  
 Reag. Range High: 1.8000  
 Blank Range Low: -0.0100  
 Blank Range High: 0.0100

Factor: -

Standard Pos: 1  
 STD-1: User Defined

STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 10-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** BUN  
**Catalog # :** B7550

Add 12ml DIH2O to 15ml sized vials. Add 40ml to 50ml sized vials.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340nm (1)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 60.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.2000  
 Point: T1  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 U/L  
 Test Range High: 80.0 U/L  
 Normal Range Low: 7.0 U/L  
 Normal Range High: 18.0 U/L

Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 2  
 Reading Last: 6  
 Reaction Limit: No  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.6000  
 Reag. Range High: 1.8000  
 Blank Range Low: -0.0100  
 Blank Range High: 0.0100

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 10-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Calcium (AZIII)  
**Catalog # :** C7529

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 600nm (5)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 320 UL  
 Start rReag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 15.0 mg/dl  
 Normal Range Low: 8.5 mg/dl  
 Normal Range High: 10.4 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: 2.000  
 Blank Range Low: No  
 Blank Range High: 2.000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Calcium (Dry)**  
**Catalog # : C7508**

Add volume DiH2O indicated on vial label. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 10 UL  
 Diluent Name: H2O  
 Volume: 10 UL

Reagent Cycle: 1  
 Volume: 360 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: 0.0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 8.5 mg/dl  
 Normal Range High: 10.4 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0500  
 Reag. Range High: 0.5000  
 Blank Range Low: -0.0800  
 Blank Range High: 0.1000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 2-03



# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Calcium**  
**Catalog # : C7503**

Prepare reagent as stated in package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 10 UL  
 Diluent Name: H2O  
 Volume: 10 UL

Reagent Cycle: 1  
 Volume: 360 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 8.5 mg/dl  
 Normal Range High: 10.4 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0500  
 Reag. Range High: 0.5000  
 Blank Range Low: -0.0800  
 Blank Range High: 0.1000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Carbon Dioxide**  
**Catalog # : C7504**

Add 10.0ml diluent to the 11ml sized vial. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mmol/l

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 2.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 mmol/l  
 Test Range High: 40 mmol/l  
 Normal Range Low: 23 mmol/l  
 Normal Range High: 34 mmol/l

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 2  
 Reading Last: 4  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -  
 Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
 recommended for calibration  
 \* USER DEFINED

It is recommended that two levels of control  
 material be assayed daily. Reorder PSI Chemistry  
 Controls Cat.# C7590-50 & C7591-50

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Chloride  
**Catalog # :** C7501

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 500nm (3)  
 Decimal Position: 0  
 Unit: mEq/L (30)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 80 mmol/l  
 Test Range High: 120 mmol/l  
 Normal Range Low: 98 mmol/l  
 Normal Range High: 106 mmol/l

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.1000  
 Reag. Range High: 0.2500  
 Blank Range Low: -0.0900  
 Blank Range High: 0.1000

Factor: -  
 Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control  
 material be assayed daily. Reorder PSI Chemistry  
 Controls Cat.# C7590-50 & C7591-50.

Rev. 6-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Cholesterol (Liquid)  
**Catalog # :** C7510

Ready to use liquid.

### GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S (1)  
Calibration Mode: Calibrator (2)  
Reagent Blank: Reag/Dil (2)  
Cleaner: No (1)

Wavelength: 500nm (3)  
Decimal Position: 0  
Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No

Sample Cycle: 2  
Volume: 3.0 UL  
Diluent Name: H2O  
Volume: 50.0 UL

Reagent Cycle: 1  
Volume: 250 UL  
Start Reag. 1 Cycle: -  
Volume: -  
Diluent: -

### CALCULATION

Sample Limit: 0.1000  
Point: T1  
Reac. Direction: Increase (1)  
Check: On

Convers. Factor: 1.00000  
Offset: 0.00000

Test Range Low: 0.0 mg/dl  
Test Range High: 500 mg/dl  
Normal Range Low: 120 mg/dl  
Normal Range High: 240 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
Reading First: 1  
Reading Last: 13  
Reaction Limit: -  
Point: -

### CALIBRATION

Calib Interval: On Request (3)  
Blank  
Reag. Range Low: -0.0700  
Reag. Range High: 0.1300  
Blank Range Low: -0.0500  
Blank Range High: 0.0500

Factor: -  
Standard Pos: 1  
STD-1: User Defined  
STD-2:  
STD-3:

Replicate: Dupl (2)  
Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
Assign:  
High:  
CS2 Pos \* Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control  
material be assayed daily. Reorder PSI Chemistry  
Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Cholesterol**  
**Catalog # : C7509**

Add 40ml distilled water to 50ml vial. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 500nm (3)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 165 UL  
 Start Reag.1 Cycle : -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.1000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 500 mg/dl  
 Normal Range Low: 120 mg/dl  
 Normal Range High: 240 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 13  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0700  
 Reag. Range High: 0.1300  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: CPK, Liquid**  
**Catalog # : C7522**

Single working reagent is prepared by mixing 4 parts R1 with 1 part R2.

Single reagent   2Part app.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)      R-S-SR1  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No  
 Sample Cycle: 1                      2  
 Volume: 5.0 UL                      5.0 UL  
 Diluent Name: H2O                      H2O  
 Volume: 30.0 UL                      15.0 UL  
 Reagent Cycle: 1                      1  
 Volume: 150 UL                      120 UL  
 Start Reag. 1 Cycle                      1  
 Volume:                                  30.0 UL  
 Diluent Name:                                  H2O  
 Diluent:                                      15.0 UL

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0500  
 Reag. Range High: 0.7000  
 Blank Range Low: -0.0900  
 Blank Range High: 0.0800  
 Factor: 9914  
 Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos      \* Low: \* (User Defined)  
 Assign:  
 High:  
 CS2 Pos      \* Low: \* (User Defined)  
 Assign:  
 High:  
 CS3 Pos:      No

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)      \*USER DEFINED  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Test Range Low: 0 U/L  
 Test Range High: 2000 U/L  
 Normal Range Low: 25 U/L  
 Normal Range High: 192 U/L

Rev. 2-03

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 5  
 Reading Last: 12  
 Reaction Limit: NO  
 Point: -

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** CPK  
**Catalog # :** C7512

Add 5ml, 12ml and 40ml to 6.5, 15 and 50ml sized vials respectfully. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1) Reag. Range Low: 0.0500  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.5000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 2000 U/L  
 Normal Range Low: 25 U/L  
 Normal Range High: 192 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 10  
 Reaction Limit: 0.235  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range High: 0.7000  
 Blank Range Low: -0.0900  
 Blank Range High: 0.0800

Factor: 8574

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** CK-MB  
**Catalog # :** C7562

Reconstitute with volume of CK-MB diluent stated on reagent vial.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340 (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 16.0 UL  
 Diluent: 10.0 UL

Reagent Cycle: 1  
 Volume: 150 UL  
 Start Reag 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 2000 U/L  
 Normal Range Low: 0 U/L  
 Normal Range High: 22 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 10  
 Reading Last: 15  
 Reaction Limit: 0.900  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0000  
 Reag. Range High: 0.7000  
 Blank Range Low: -0.0020  
 Blank Range High: 0.0800

Factor: 6200

Calibrator Cup Pos.-  
 CAL-1: -

Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\*USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI CK-MB Controls Cat.# C7562-CTL

Rev. 2-03



# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Creatinine CS  
 Catalog # : C7539

### 2-Part reagent application

#### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Slope Avg.  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before

Wavelength: 500nm (3)  
 Decimal Position: 1  
 Unit: mg/dl (12)

#### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 15.0 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 100 UL  
 Start Reag. 1 Cycle: 2  
 Volume: 50.0 UL  
 Diluent: H2O  
 Volume: 10.0 UL

#### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 25.0 mg/dl  
 Normal Range Low: 0.4 mg/dl  
 Normal Range High: 1.4 mg/dl

Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 3  
 Reading Last: 7  
 Reaction Limit: No (1)  
 Point: -

#### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0100  
 Reag. Range High: 0.45000  
 Blank Range Low: -0.0030  
 Blank Range High: 0.0300

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

#### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
 recommended for calibration.

\* USER DEFINED.

It is recommended that two levels of control  
 material be assayed daily. Reorder PSI Chemistry  
 Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** CRP (HS)  
**Catalog # :** C7564

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR (1)  
 Calibration Mode: STD NonLinear  
 Reagent Blank: No  
 Cleaner: No

Wavelength: 550 nm  
 Decimal Position: 2  
 Unit: mg/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 3  
 Volume: 15.0 UL

Diluent Name: Diluent  
 Volume: 30.0 UL (use rgt 1)

Reagent Cycle: 1  
 Volume: 130 UL  
 Start Reag 1 Cycle: 1  
 Volume: 90.0 U/L  
 Diluent Name: use rgt 1  
 Diluent: 5.0U/L

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: No  
 Test Range High: No  
 Normal Range Low: \*  
 Normal Range High: \*

Number of Steps: 1

Calc. Step A: Kinetic  
 Reading First: 3  
 Reading Last: 10  
 Reaction Limit: No  
 Point: -

### CALIBRATION

Calib Interval: On Request  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -  
 Standard Pos. \*  
 Std-1: User Defined  
 Std-2:  
 Std-3:

Replicate: Single  
 Deviation: No  
 Calc. Model: Logit/Logit4  
 Correction Std: No

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* User defined.  
 Use saline as 0.0 standard  
 Enter calibrator values.

Rev. 11-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Fructosamine**  
**Catalog # : F7546**

Prepare reagent according to package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope Avg (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 550nm (4)  
 Decimal Position: 0  
 Unit: mmol/L

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 10 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 200 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 1.2  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On (1)

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mmol/L  
 Test Range High: 8.5 mmol/L  
 Normal Range Low: 1.3 mmol/L  
 Normal Range High: 2.85 mmol/L

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 20  
 Reading Last: 33  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: 0.15  
 Blank Range Low: -0.003  
 Blank Range High: 0.003

Factor: -

Standard Pos: -  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 5%

### Control

CS1 Pos \* Low: User Defined  
 Assign:  
 High:  
 CS2 Pos \* Low: User Defined  
 Assign:  
 High:  
 CS3 Pos: NO

(Non-validated application)

\* USER DEFINED

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Gamma GT (Liquid)  
**Catalog # :** G7571

Single working reagent is prepared by mixing 4 parts R1 with 1 part R2.

	<u>Single Reagent</u>	<u>2Part app.</u>		
<b><u>GENERAL</u></b>			<b><u>CALIBRATION</u></b>	
Measurement Mode:	Absorb		Calib Interval:	On Request (3)
Reaction Mode:	R-S (1)	R-S-SR1	Blank	
Calibration Mode:	Factor (1)		Reag. Range Low:	0.1000
Reagent Blank:	Reag/Dil (2)		Reag. Range High:	0.8000
Cleaner:	No (1)		Blank Range Low:	-0.0100
			Blank Range High:	0.0100
Wavelength:	405nm (2)		Factor:	3938
Decimal Position:	0		Standard Pos:	-
Unit:	U/L (21)		STD-1:	-
			STD-2:	-
<b><u>ANALYSIS</u></b>			STD-3:	-
Post Dil. Factor:	No		Replicate:	-
Post Conc. Factor:	No		Deviation:	-
Sample Cycle:	2	2		
Volume:	10.0 UL	10.0 UL	<b><u>Control</u></b>	
Diluent Name:	H2O	H2O	CS1 Pos	* Low: (User Defined)
Volume:	20.0 UL	10 UL		Assign:
Reagent Cycle:	1	1		High:
Volume:	165 UL	130 UL	CS2 Pos	* Low: (User Defined)
Start Reag 1 Cycle:	-	1		Assign:
Volume:	-	35 UL		High:
Diluent Name:	-	H2O	CS3 Pos:	No
Volume:	-	10 UL		
<b><u>CALCULATION</u></b>				
Sample Limit:	NO			* USER DEFINED
Point:	-			
Reac. Direction:	Increase (1)			It is recommended that two levels of control
Check:	On			material be assayed daily. Reorder PSI Chemistry
Convers. Factor:	1.00000			Controls Cat.# C7590-50 & C7591-50.
Offset:	0.00000			Rev. 2-03
Test Range Low:	0.0 U/L			
Test Range High:	1000 U/L			
Normal Range Low:	8 U/L			
Normal Range High:	54 U/L			
Number of Steps:	1			
Calc. Step A:	Kinsearch (3)			
Reading First:	3			
Reading Last:	11			
Reaction Limit:	NO			
Point:	-			

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Gamma GT Soluble  
**Catalog # :** G7570

Add 10.0ml DH2O to 10ml sized vial.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 405nm (2)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 10.0 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 165 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.5000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 1000 mg/dl  
 Normal Range Low: 8 mg/dl  
 Normal Range High: 54 mg/dl

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 11  
 Reaction Limit: 0.300  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.1000  
 Reag. Range High: 0.8000  
 Blank Range Low: -0.0100  
 Blank Range High: 0.0100

Factor: 3282

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Glucose Hex (Liquid)  
**Catalog # :** G7517

Reagent supplied as ready to use liquid.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 600 mg/dl  
 Normal Range Low: 65 mg/dl  
 Normal Range High: 110 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: CB  
 Reading Last: 13  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.050  
 Reag. Range High: 0.3000  
 Blank Range Low: -0.050  
 Blank Range High: 0.3000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 10-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Glucose Hex  
**Catalog # :** G7518

Add 12ml and 40.0ml to 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 200 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.3000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 600 mg/dl  
 Normal Range Low: 65 mg/dl  
 Normal Range High: 110 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 6  
 Reaction Limit: -  
 Point:

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0030  
 Reag. Range High: 0.3000  
 Blank Range Low: -0.0700  
 Blank Range High: 0.0600

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: HDL Cholesterol  
 Catalog # : H7507/H7511

See package insert.

NOTE: Cholesterol Reagent should be reconstituted as per cholesterol Mira application instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 500nm (3)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 6.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 165 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.2000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 250 mg/dl  
 Normal Range Low: 30 mg/dl  
 Normal Range High: 75 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 13  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0700  
 Reag. Range High: 0.1300  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -

Standard Pos: 2  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

NOTE: Final results must be multiplied by 2 or the calibrator value entered into the instrument doubled.  
 Cholesterol calibrator catalog# C7574-50 is recommended for calibration.

\* USER DEFINED

Rev. 2-03



# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Hemoglobin**  
**Catalog # : H7504**

Ready to use liquid.

**GENERAL**

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 550nm  
 Decimal Position: 1  
 Unit: g/dl

**ANALYSIS**

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

**CALCULATION**

Sample Limit: 0.3000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 g/dl  
 Test Range High: 25.0 g/dl  
 Normal Range Low: 11.0 g/dl  
 Normal Range High: 25.0 g/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 13  
 Reaction Limit: -  
 Point: -

**CALIBRATION**

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0700  
 Reag. Range High: 0.1300  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

**Control**

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Hemoglobin Standard catalog #H7504-STD recommended for calibration.

\*USER DEFINED

\*\*Non-validated application

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Total Iron  
**Catalog # :** I7504 / I7505

Reagents are supplied ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1                      Blank  
 Calibration Mode: Std Lin or Slope Avg.  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 0  
 Unit: ug/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 50 UL  
 Diluent Name: H2O  
 Volume: 40 UL

Reagent Cycle: 1  
 Volume: 240 UL\*  
 Start Reag 1 Cycle: 4  
 Volume: 5 UL\*\*  
 Diluent: 0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 ug/dl  
 Test Range High: 510 ug/dl  
 Normal Range Low: 60 ug/dl  
 Normal Range High: 150 ug/dl

Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 3  
 Reading Last: 28  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Reag. Range Low: No  
 Reag. Range High: 2.000  
 Blank Range Low: No  
 Blank Range High: 2.000

Factor: -

Standard Pos:  
                   STD-1: User Defined  
                   STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos           # Low: (User Defined)  
                           Assign:  
                           High:  
 CS2 Pos           # Low: (User Defined)  
                           Assign:  
                           High:  
 CS3 Pos:           No

\*Iron Buffer Reagent  
 \*\* Iron Color Reagent

### # USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistr Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Lactate  
**Catalog # :** L7596

Prepare reagent according to package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope Avg. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 600nm (5)  
 Decimal Position: 1  
 Unit: mmol/L

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Rec. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 mmol/L  
 Test Range High: 20 mmol/L  
 Normal Range Low: 0.5 mmol/L  
 Normal Range High: 2.2 mmol/L

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 13  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.100  
 Reag. Range High: 2.000  
 Blank Range Low: -0.1000  
 Blank Range High: 1.5000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** HBDH  
**Catalog # :** H7569

Add 5ml, 12ml to 6.5, 15 ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 150 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: NO  
 Point: T1  
 Reac. Direction: Decrease (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 800 U/L  
 Normal Range Low: 68 U/L  
 Normal Range High: 135 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 5  
 Reading Last: 12  
 Reaction Limit: No  
 Point:

### CALIBRATION

Calib Interval: On Request (3)  
 Blank

Reag. Range Low: 0.4000  
 Reag. Range High: 1.9000  
 Blank Range Low: -0.0500  
 Blank Range High: 0.5000

Factor: 9914

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

Control

CS1 Pos \* Low: (User Defined)  
 Assign: High:

CS2 Pos \* Low: (User Defined)  
 Assign: High:

CS3 Pos: No

\* USER DEFINED

\*\*non-validated application

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** LDH  
**Catalog # :** L7535

Add 5ml, 12ml and 40ml to 6.5, 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Select (3) After AST,ALT

Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.5000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 1250 U/L  
 Normal Range Low: 80 U/L  
 Normal Range High: 285 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 11  
 Reaction Limit: 0.1500  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.1000  
 Reag. Range High: 0.7000  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: 8574

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

Control  
 CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

### \* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** LDH, Liquid  
**Catalog # :** L7572

Prepare single working reagent by mixing 4 parts R1 and 1 part R2.

Single Reagent      2Part app.

**GENERAL**

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)      R-S-SR1  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

**ANALYSIS**

Post Dil. Factor: No  
 Post Conc. Factor: No  
 Sample Cycle: 2                      2  
 Volume: 5.0 UL                      5.0 U/L  
 Diluent Name: H2O                      H2O  
 Volume: 30.0 UL                      15.0 UL  
 Reagent Cycle: 1                      1  
 Volume: 150UL                      120 UL  
 Start Reag 1 Cycle: -                      1  
 Volume: -                      30 UL  
 Diluent Name: -                      H2O  
 Diluent: -                      15 UL

**CALCULATION**

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On  
 Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 1400 U/L  
 Normal Range Low: 80 U/L  
 Normal Range High: 285 U/L

Number of Steps: 1

Calc. Step A: Kinsearch (3)  
 Reading First: 5  
 Reading Last: 12  
 Reaction Limit: No  
 Point: -

**CALIBRATION**

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.1000  
 Reag. Range High: 0.4000  
 Blank Range Low: -0.0050  
 Blank Range High: 0.0010  
 Factor: 9914  
 Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Replicate: -  
 Deviation: -

**Control**

CS1 Pos                      \*Low: (User Defined)  
    Assign:  
    High:  
 CS2 Pos                      \*Low: (User Defined)  
    Assign:  
    High:  
 CS3 Pos:                      No

\*USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03

# ROCHE COBAS MIRA

## LD-P (LDH) – CATALOG# L7536

Add 5ml, 12ml and 40ml to 6, 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Select (3) After AST,ALT  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No  
 Sample Cycle: 1  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL  
 Reagent Cycle: 1  
 Volume: 125 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.5000  
 Point: T1  
 Reac. Direction: Decrease  
 Check: On  
 Convers. Factor: 1.00000  
 Offset: 0.00000  
 Test Range Low: 0 U/L  
 Test Range High: 1200 U/L  
 Normal Range Low: 190 U/L  
 Normal Range High: 560 U/L  
 Number of Steps: 1  
 Calc. Step A: Kinsearch (3)  
 Reading First: 3  
 Reading Last: 11  
 Reaction Limit: 0.1500  
 Point: T1

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.500  
 Reag. Range High: 1.900  
 Blank Range Low: -0.500  
 Blank Range High: 0.500  
 Factor: 8574  
 Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:  
 Replicate: -  
 Deviation: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

Rev: 5/01

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** Lipase  
**Catalog # :** L7503

Add 10ml substrate buffer to lipase substrate. Use activator reagent as start reagent.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 0  
 Unit: U/L (21)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 3.0 UL  
 Diluent: 10.0 UL

Reagent Cycle: 1  
 Volume: 180 UL  
 Start Reag 1 Cycle: 8  
 Volume: 60.0 UL  
 Diluent: 10.0 UL

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 600 U/L  
 Normal Range Low: 0 U/L  
 Normal Range High: 62 U/L  
 Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 14  
 Reading Last: 20  
 Reaction Limit: No (1)  
 Point: -

### CALIBRATION

Calib Interval: On Request  
 Blank  
 Reag. Range Low: -0.1000  
 Reag. Range High: 0.5000  
 Blank Range Low: -0.0900  
 Blank Range High: 0.1000

Factor: -  
 Calibrator Pos: 3  
 CAL-1: User Defined  
 CAL-2:  
 CAL-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Use Lipase Standard as Calibrator.

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev. 2-03



# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Magnesium  
 Catalog # : M7527

Reagent is supplied ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator. (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: mg/dl

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 15.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 6.0 mg/dl  
 Normal Range Low: 1.6 mg/dl  
 Normal Range High: 3.0 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 4  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.1000  
 Reag. Range High: 0.7000  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -  
 Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Chemistry Calibrator catalog #C7506-50 recommended for calibration

Rev: 7/02

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: Microalbumin (2-Point)  
Catalog # : M7562

Reagents are supplied ready to use.

## GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S-SR1  
Calibration Mode: Slope Avg. (2)  
Reagent Blank: No  
Cleaner: No  
  
Wavelength: 340nm  
Decimal Position: 2  
Unit: mg/dl (12)

## ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No  
  
Sample Cycle: 1  
Volume: 12.0 UL  
Diluent Name: H2O  
Volume: 10.0 UL  
  
Reagent Cycle: 1  
Volume: 270 UL  
Start Reag 1 Cycle: 8  
Volume: 90.0 UL  
Diluent: 0.0

## CALCULATION

Sample Limit: No  
Point: -  
Reac. Direction: Increase (1)  
Check: Off

Convers. Factor: 1.00000  
Offset: 0.00000

Test Range Low: No  
Test Range High: No  
Normal Range Low: No  
Normal Range High: No

Number of Steps: 1  
Calc. Step A: Endpoint (1)  
Reading First: 7  
Reading Last: 22  
Reaction Limit: -  
Point: -

## CALIBRATION

Calib Interval: On Request (3)  
Blank  
Reag. Range Low: No  
Reag. Range High: No  
Blank Range Low: No  
Blank Range High: No  
  
Factor: -  
  
Standard Pos: 2  
STD-1: User Defined  
STD-2:  
STD-3:

Replicate: Single  
Deviation: No

## Control

CS1 Pos \* Low: (User Defined)  
Assign:  
High:  
CS2 Pos \* Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

Use saline as 0.0 standard  
Microalbumin calibrator set catalog# M7562-CAL  
required for calibration.

\* USER DEFINED

Rev: 11-03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Microalbumin (multi-Point)  
 Catalog # : M7562

Reagents are supplied ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Lin Inter  
 Reagent Blank: No  
 Cleaner: No

Wavelength: 340nm  
 Decimal Position: 2  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 12.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 270 UL  
 Start Reag 1 Cycle: 8  
 Volume: 90.0 UL  
 Diluent: 0.0

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off  
 Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: No  
 Test Range High: No  
 Normal Range Low: No  
 Normal Range High: No

Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 7  
 Reading Last: 22  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -

Standard Pos: 2  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Single  
 Deviation: No

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

Use saline as 0.0 standard.  
 Microalbumin calibrator set catalog# M7562-CAL  
 required for calibration.

\* USER DEFINED

Rev: 11-03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Microprotein  
 Catalog # : P7582

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 600nm (5)  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 4 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 200 mg/dl  
 Normal Range Low: \*  
 Normal Range High: \*

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: CB  
 Reading Last: 24  
 Reaction Limit: No  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -

Standard Pos: \*  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Triplicate (3)  
 Deviation: No

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\*User Defined

Rev: 2-03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Phosphorus  
 Catalog # : P7516

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 340nm (1)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 7.0 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 300 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 mg/dl  
 Test Range High: 15.0 mg/dl

Normal Range Low: 2.5 mg/dl  
 Normal Range High: 4.8 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: CB  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.1500  
 Reag. Range High: 0.3900  
 Blank Range Low: -0.0500  
 Blank Range High: 0.1000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
 recommended for calibration.

\* USER DEFINED  
 Rev: 2-03

It is recommended that two levels of control  
 material be assayed daily. Reorder PSI Chemistry  
 Controls Cat.# C7590-50 & C7591-50.

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Total Protein  
 Catalog # : T7528

Reagents are ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 550nm (4)  
 Decimal Position: 1  
 Unit: g/dl (11)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 5.0 UL  
 Diluent Name: H2O  
 Volume: 20.0 UL

Reagent Cycle: 1  
 Volume: 250 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.3000  
 Point: 1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0.0 g/dl  
 Test Range High: 15.0 g/dl  
 Normal Range Low: 6.2 g/dl  
 Normal Range High: 8.5 g/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 5  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.1000  
 Reag. Range High: 0.4000  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

\* USER DEFINED

Rev: 2-03

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: Rheumatoid Factor  
Catalog # : R7568

Reagents are supplied ready to use.

## GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S-SR1  
Calibration Mode: Lin Inter  
Reagent Blank: No  
Cleaner: No

Wavelength: 340nm  
Decimal Position: 0  
Unit: U/ML

## ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No

Sample Cycle: 1  
Volume: 15.0 UL  
Diluent Name: H2O  
Volume: 10.0 UL

Reagent Cycle: 1  
Volume: 250 UL  
Start Reag 1 Cycle: 12  
Volume: 50.0 UL  
Diluent: 10.0 UL

## CALCULATION

Sample Limit: No  
Point: -  
Reac. Direction: Increase (1)  
Check: Off

Convers. Factor: 1.00000  
Offset: 0.00000

Test Range Low: No  
Test Range High: No  
Normal Range Low: No  
Normal Range High: No

Number of Steps: 1  
Calc. Step A: Endpoint (1)  
Reading First: 11  
Reading Last: 22  
Reaction Limit: -  
Point: -

## CALIBRATION

Calib Interval: On Request (3)

Blank  
Reag. Range Low: No  
Reag. Range High: No  
Blank Range Low: No  
Blank Range High: No

Factor: -

### Standard Pos:

STD-1: User Defined  
STD-2: 2-6 USER DEFINED  
STD-3:

Replicate: Single  
Deviation: No

## Control

CS1 Pos \* Low: (User Defined)  
Assign:  
High:  
CS2 Pos \* Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

Use saline as 0.0 standard.  
RF calibrator set catalog# R7568-CAL required for calibration.

\* USER DEFINED

Rev: 11/03

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Sodium  
 Catalog # : S7571

Prepare reagent according to package insert instructions.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: No (1)

Wavelength: 550nm (4)  
 Decimal Position: 0  
 Unit: mmol/L

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 10 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Reagent Cycle: 1  
 Volume: 250 UL  
 Start Reag. 1 Cycle: 4  
 Volume: 95 UL  
 Diluent:

### CALCULATION

Sample Limit: No  
 Point: -  
 Reac. Direction: Increase (1)  
 Check: Off

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: No  
 Test Range High: No  
 Normal Range Low: No  
 Normal Range High: No  
 Number of Steps: 1  
 Calc. Step A: Endpoint (1)  
 Reading First: 3  
 Reading Last: 7  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: Each Run  
 Blank  
 Reag. Range Low: No  
 Reag. Range High: No  
 Blank Range Low: No  
 Blank Range High: No

Factor: -

Standard Pos: 7  
 STD-1: Standard  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: No

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:

CS3 Pos: No

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

\* USER DEFINED  
 Rev: 5/01

Serum samples should be prepared as stated in package insert. R1 is the Acid Reagent. R2 is the Color Reagent.



# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: Triglyceride-GPO (Liq)  
 Catalog # : T7532

Reagent provided ready to use.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 500nm  
 Decimal Position: 0  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 2  
 Volume: 3.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 200 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.6000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 mg/dl  
 Test Range High: 1000 mg/dl  
 Normal Range Low: 36 mg/dl  
 Normal Range High: 165 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: 1  
 Reading Last: 7  
 Reaction Limit: -  
 Point:

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0700  
 Reag. Range High: 0.3000  
 Blank Range Low: -0.0900  
 Blank Range High: 0.0700

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
 recommended for calibration.

It is recommended that two levels of control  
 material be assayed daily. Reorder PSI Chemistry  
 Controls Cat.# C7590-50 & C7591-50.

\* USER DEFINED  
 Rev: 2-03

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: Triglyceride-GPO  
Catalog # : T7531

Add 12ml and 40ml to 15ml and 50ml vials respectively. Swirl to dissolve.

## GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S (1)  
Calibration Mode: Calibrator (2)  
Reagent Blank: Reag/Dil (2)  
Cleaner: Before (2)

Wavelength: 550nm  
Decimal Position: 0  
Unit: mg/dl (12)

## ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No

Sample Cycle: 1  
Volume: 3.0 UL  
Diluent Name: H2O  
Volume: 50.0 UL

Reagent Cycle: 1  
Volume: 200 UL  
Start Reag. 1 Cycle: -  
Volume: -  
Diluent: -

## CALCULATION

Sample Limit: 0.6000  
Point: T1  
Reac. Direction: Increase (1)  
Check: On

Convers. Factor: 1.00000  
Offset: 0.00000

Test Range Low: 0 mg/dl  
Test Range High: 1000 mg/dl  
Normal Range Low: 36 mg/dl  
Normal Range High: 165 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
Reading First: T1  
Reading Last: 7  
Reaction Limit: -  
Point: -

## CALIBRATION

Calib Interval: On Request (3)  
Blank  
Reag. Range Low: -0.0700  
Reag. Range High: 0.3000  
Blank Range Low: -0.0900  
Blank Range High: 0.0700

Factor: -  
Standard Pos: 1  
STD-1: User Defined  
STD-2:  
STD-3:

Replicate: Dupl (2)  
Deviation: 10%

## Control

CS1 Pos \* Low: (User Defined)  
Assign:  
High:  
CS2 Pos \* Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

Chemistry Calibrator catalog #C7506-50  
recommended for calibration.

It is recommended that two levels of control  
material be assayed daily. Reorder PSI Chemistry  
Controls Cat.# C7590-50 & C7591-50.

\* USER DEFINED  
Rev: 2-03

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: UIBC - Notes  
Catalog # : I7504, I7506

**Calculations:**

Total Iron Binding Capacity is calculated as follows:

$$\text{TIBC (ug/dl)} = \text{Total Iron} + \text{UIBC}$$

% Saturation is calculated as Follows:

$$\% \text{ Saturation} = \frac{\text{Serum Iron} \times 100}{\text{TIBC}}$$

**Linearity:**

Up to 500 ug/dl

If UIBC in the sample exceeds 500ug/dl, dilute 1 part sample with 2 parts saline and reassay. Multiply result by 3 to compensate for dilution.

**Procedure:**

1. Fill reagent containers with sufficient amounts of reagents.
2. Enter assay parameters as outlined.
3. Assign the 500 standard to a calibrator position as outline in the instrument settings and place the standard in that position. Calibrate the assay by requesting a pre-calibration (PC) to establish the appropriate calibration factor.
4. After a calibration factor has been determined modify the instrument settings as follows.

**GENERAL**

Calibration Mode : Factor.....1

**CALCULATION**

Reaction Direction: Decrease.....2

**CALIBRATION**

Factor: .....(Enter Calibration Factor from step 3)

5. Start Assay.

**Procedural Notes:**

1. The cup position and calibrator concentration will be automatically deleted when the calibrator mode is changed to Factor. The calibration interval is set to each run so that the reagent blank will be updated with each run. The reagent blank will be updated automatically by the instrument using the instrument diluent as sample. Calibrators are not required to update the blank.
2. It is not necessary to establish a calibration factor with each run. However, the calibration factor should be verified with each reagent lot change, instrument maintenance, or if control material indicates the need for recalibration.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: UIBC  
Catalog # : I7504, I7506

## GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S-SR1-SR2  
Calibration Mode: Slope Avg  
Reagent Blank: Reag/Dil (2)  
Cleaner: Before

Wavelength: 550nm (4)  
Decimal Position: 0  
Unit: ug/dl

## ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No

Sample Cycle: 1  
Volume: 25 UL  
Diluent Name: H2O  
Volume: 10 UL

Reagent Cycle: 1  
Volume: 170 UL  
Start Reag 1 Cycle: 2  
Volume: 30 UL  
Diluent: 10 UL  
Start Reag 2 Cycle: 10  
Volume: 5 UL  
Diluent: 10

## CALCULATION

Sample Limit: No  
Point: -  
Reac. Direction: Increase  
Check: On

Convers. Factor: 1.00000  
Offset: 0.00000

Test Range Low: 0 ug/dl  
Test Range High: 510 ug/dl  
Normal Range Low: 130 ug/dl  
Normal Range High: 375 ug/dl

Number of Steps: 1  
Calc. Step A: Endpoint (1)  
Reading First: 9  
Reading Last: 24  
Reaction Limit: -  
Point: -

## CALIBRATION

Calib Interval: Each Run  
Blank  
Reag. Range Low: -0.010  
Reag. Range High: 2.000  
Blank Range Low: -0.010  
Blank Range High: 2.000

Factor: enter Iron calibration factor

Standard Pos: User Defined  
STD-1: 500  
STD-2:  
STD-3:

Replicate: Dupl (2)  
Deviation: 10%

Control  
CS1 Pos # Low: (User Defined)  
Assign:  
High:  
CS2 Pos # Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

**Reagent 1:** UIBC buffer **Start reagent 1:** Iron Standard  
**Start reagent 2:** Iron Color

# USER DEFINED

See Notes for detailed instructions on use of reagent and application

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Rev: 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Uric Acid (Liquid)**  
**Catalog # : U7581**

Reagent supplied as a ready to use product.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Slope AVG (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: no  
  
 Wavelength: 500nm (3)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 4.0 UL  
 Diluent Name: H2O  
 Volume: 30.0 UL

Reagent Cycle: 1  
 Volume: 150 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: no  
 Point:  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 2.2 mg/dl  
 Normal Range High: 7.7 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: CB  
 Reading Last: 15  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: 0.0000  
 Reag. Range High: 0.2500  
 Blank Range Low: -0.0500  
 Blank Range High: 0.1000

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

\* USER DEFINED

Rev: 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer: Cobas Mira**  
**Test: Uric Acid**  
**Catalog # : U7580**

Add 12ml and 40ml to 15 and 50ml sized vials respectively. Swirl to dissolve.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Calibrator (2)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)

Wavelength: 500nm (3)  
 Decimal Position: 1  
 Unit: mg/dl (12)

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 6.0 UL  
 Diluent Name: H2O  
 Volume: 50.0 UL

Reagent Cycle: 1  
 Volume: 200 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### CALCULATION

Sample Limit: 0.6000  
 Point: T1  
 Reac. Direction: Increase (1)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 mg/dl  
 Test Range High: 20.0 mg/dl  
 Normal Range Low: 2.2 mg/dl  
 Normal Range High: 7.7 mg/dl

Number of Steps: 1

Calc. Step A: Endpoint (1)  
 Reading First: T1  
 Reading Last: 8  
 Reaction Limit: -  
 Point: - Rev: 5/01

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.0700  
 Reag. Range High: 0.3500  
 Blank Range Low: -0.0500  
 Blank Range High: 0.0500

Factor: -

Standard Pos: 1  
 STD-1: User Defined  
 STD-2:  
 STD-3:

Replicate: Dupl (2)  
 Deviation: 10%

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

Chemistry Calibrator catalog #C7506-50 recommended for calibration

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

\* USER DEFINED

Rev: 2-03

# Pointe Scientific, Inc.

## Instrument Application

**Analyzer:** Cobas Mira  
**Test:** G6PD  
**Catalog # :** G7583

**Reagent preparation:** Prepare working reagent by adding 6 mls DH2O to the stated 6 ml vial and let dissolve. Now add 12 mls of the R2 reagent to the same vial. This is your working reagent.

**Sample preparation:** Add 100 ul whole blood to 0.9 mls lyse reagent and let stand 5 minutes. Mix well.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S (1)  
 Calibration Mode: Factor (1)  
 Reagent Blank: Reag/Dil (2)  
 Cleaner: Before (2)  
 Wavelength: 340nm (1)  
 Decimal Position: 0  
 Unit: U/L (21)

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.010  
 Reag. Range High: 0.240  
 Blank Range Low: -0.010  
 Blank Range High: 0.010  
 Factor: 98377

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Standard Pos: -  
 STD-1: -  
 STD-2:  
 STD-3:

Sample Cycle: 1  
 Volume: 7.0 UL  
 Diluent Name: H2O  
 Volume: 10.0 UL

Replicate: -  
 Deviation: -

Reagent Cycle: 1  
 Volume: 240 UL  
 Start Reag. 1 Cycle: -  
 Volume: -  
 Diluent: -

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

### CALCULATION

Sample Limit: NO  
 Point: T1  
 Reac. Direction: Increase (2)  
 Check: On

Calculations:  $G6PDH (U/gHgb) = G6PDH (U/l) / (10 \times Hgb)$

\* USER DEFINED  
 Rev: 2-03

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0 U/L  
 Test Range High: 3000 U/L  
 Normal Range Low: NO  
 Normal Range High: NO

Number of Steps: 1

Calc. Step A: Kinetic (2)  
 Reading First: 8  
 Reading Last: 19  
 Reaction Limit: .1000  
 Point: T1

# Pointe Scientific, Inc.

## Instrument Application

Analyzer: Cobas Mira  
 Test: HbA1c  
 Catalog # : H7541

**Reagent preparation:** Reagent provided ready to use.

**Sample preparation:** Add 10 ul packed cells to 1.0 mls hemolysis reagent and let stand 5 minutes. Mix well.

### GENERAL

Measurement Mode: Absorb  
 Reaction Mode: R-S-SR1  
 Calibration Mode: Std. NONLIN  
 Reagent Blank: NO Blk  
 Cleaner: NO  
 Wavelength: 600nm (1)  
 Decimal Position: 2  
 Unit: %

### ANALYSIS

Post Dil. Factor: No  
 Post Conc. Factor: No

Sample Cycle: 1  
 Volume: 7.0 UL  
 Diluent Name: H2O  
 Volume: 3.0 UL

Reagent Cycle: 1  
 Volume: 240 UL  
 Start Reag. 1 Cycle: 12  
 Volume: 80 UL  
 Diluent: 5.0 UL

### CALCULATION

Sample Limit: NO  
 Point: -  
 Reac. Direction: Increase (2)  
 Check: On

Convers. Factor: 1.00000  
 Offset: 0.00000

Test Range Low: 0  
 Test Range High: 16.0  
 Normal Range Low: NO  
 Normal Range High: NO

Number of Steps: 1

Calc. Step A: ENDPOINT  
 Reading First: 13  
 Reading Last: 24  
 Reaction Limit: -  
 Point: -

### CALIBRATION

Calib Interval: On Request (3)  
 Blank  
 Reag. Range Low: -0.010  
 Reag. Range High: 0.240  
 Blank Range Low: -0.010  
 Blank Range High: 0.010  
 Factor: -

Standard Pos: 1  
 STD-1: \* STD-4: \*  
 STD-2: \* STD-5: \*  
 STD-3: \*

Replicate: Single  
 Deviation: NO  
 CALC MODEL: Logit/LOG

### Control

CS1 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS2 Pos \* Low: (User Defined)  
 Assign:  
 High:  
 CS3 Pos: No

\* USER DEFINED  
 Rev: 11/03



# Pointe Scientific, Inc.

Instrument Application

Analyzer: Cobas Mira  
Test: Lp(a)  
Catalog # : L7597

**Reagent preparation:** Reagents provided as ready to use liquids.

## GENERAL

Measurement Mode: Absorb  
Reaction Mode: R-S -SR1  
Calibration Mode: Std. Nonlin  
Reagent Blank: Reag/Dil (2)  
Cleaner: NO  
Wavelength: 340nm  
Decimal Position: 1  
Unit: mg/dL

## ANALYSIS

Post Dil. Factor: No  
Post Conc. Factor: No  
  
Sample Cycle: 1  
Volume: 15.0 uL  
Diluent Name: H2O  
Volume: 20.0 uL  
  
Reagent Cycle: 1  
Volume: 300 uL  
Start Reag. 1 Cycle: 13  
Volume: 50 uL  
Diluent: H2O  
Volume: 10.0 uL

## CALCULATION

Sample Limit: NO  
Point: -  
Reac. Direction: Increase (2)  
Check: Off  
  
Convers. Factor: 1.00000  
Offset: 0.00000  
  
Test Range Low: No  
Test Range High: No  
Normal Range Low: 0.02  
Normal Range High: 0.27

Number of Steps: 1  
  
Calc. Step A: Endpoint  
Reading First: 12  
Reading Last: 25  
Reaction Limit: NO  
Point: -

## CALIBRATION

Calib Interval: On Request (3)  
Blank  
Reag. Range Low: NO  
Reag. Range High: NO  
Blank Range Low: NO  
Blank Range High: NO  
Factor: -  
  
Standard Pos: \*1  
STD-1: \*0.0  
STD-2: \*2  
STD-3: \*2

Replicate: Triple  
Deviation: -20%

## Control

CS1 Pos \* Low: (User Defined)  
Assign:  
High:  
CS2 Pos \* Low: (User Defined)  
Assign:  
High:  
CS3 Pos: No

\*1 Position of Calibrator  
\*2 Assigned Value of Calibrator

Rev. 3/03